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## 1. Warnings and Safety Instructions

### **Please note:**

*Please read the user manual and safety instructions carefully. Each person working in the laboratory must be familiar with these safety instructions prior to using the device. The manual must be accessible at any time. Not every possible safety issue is covered in this manual. Each applicant of the device is solely responsible to maintain safety and health standards and to determine and follow restrictions prior to using the device.*

*- All general safety instructions of the laboratory such as wearing protective clothing, safety glasses and gloves appropriate for the application of certain fluids must be followed at all times.*

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- Instructions provided by the reagent manufacturer must be observed at any time.

- This device is only intended for dispensing/ titrating fluids. Certain applications are either prohibited or restricted. These instructions must be strictly adhered to. (If in doubt contact the manufacturer).

- Prior to use, all connections must be checked for leakage. The correct positioning of the discharge spout and the suction tube must be verified. If the discharge spout is not secured properly fluid could leak or splatter out during the operation of the device.

- The closure cap must be removed from the discharge tube prior to dispensing/titrating. Do not start operation as long as the closure cap is still positioned on the device.

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- Make sure that you do not endanger yourself or any other person. Never work towards another person. Avoid splattering. Use the appropriate receiving vessels. Always dispense/titrate onto the inner wall of the receiving container instead of its content or onto the bottom of the container.

- Never use force!

- As soon as a dispensing/ titration has been completed position closure cap on discharge tube. Caution! Reagent may leak out!

- Use original accessories and replacement parts only.

- Once sterilization has been completed (refer to section 12) let the device cool down to room temperature before use.

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- In the event of errors, such as
  - slow moving or jammed piston
  - sticky or leaking suction/discharge recirculation valve
  - leakage
  - broken parts

stop working immediately. Before you continue to use the device read the instructions for cleaning of the device (section 11) and troubleshooting (section 17). If necessary send the unit to the manufacturers' service center for repair. Please read section 15 for more information.

- Do not place the device or the charger unit into water or other fluids, in order to avoid the risk of an electric shock!
- Make sure that the voltage of the power supply corresponds with the indicated voltage on the charger unit.

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*- Make sure that a charger unit is not used in the event that the cable or the charger itself is damaged.*

*- Never touch the charger unit's plug with wet hands!*

*- Protect the charger unit from humidity! Only use it indoors.*

### **Warning!**


*Opening the device results in loss of liability for the device and any damage caused by it.*

*- Avoid improper use.*

### **Warning!**

*To prevent damage to the instrument, always use manufacturer's original charger unit, accessories and spares.*

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### 3. Functional principle

- Combined with the appropriate electronic module, the instrument is a battery-operated motorized system for dispensing or titrating fluids. The basic instrument (with thread A45) with its variety of adapters (A32, A38, S40) can be mounted on common reagent bottles.

- The system operates on the piston stroke principle with check valves. The piston unit is belt-driven by a DC motor with incremental step control via micro-processor. Power is supplied by 8 AA size batteries.

#### **Note:**

The upper part of the instrument is not impermeable to gas and humidity. Therefore it must not be used in aggressive vapors.

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- When the piston unit is raised, the suction valve opens. From a reservoir bottle, fluid is drawn into the cylinder through the suction tube. When the suction process is finished, the suction valve closes.

- When the piston unit moves down, the discharge valve opens and releases the fluid in the cylinder.

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#### **4. Prohibited applications**

*The device may not be used with:*

- *Fluids that attack Halar (ECTFE), FEP and Hastelloy*
- *Suspensions, since solid particles may block the valves*
- *Hydrofluoric acid, since it attacks glass*
- *Easily flammable fluids (flammable lower than 0°C, e.g. ether, acetone)*
- *Use in strictly Ex-protected areas (Ex-area zone 1)*

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## **5. Restricted applications**

- The range of the working temperature for the device and the fluid is in between +4°C and +40°C.
- Use of the device only in areas where a dangerous and explosive atmosphere can be ruled out.
- Charge only in Ex-secured areas.
- Do not place unattended into Ex-protected areas.

### **Important! Please note!**

Never use force to release a jammed piston. Immediately stop operation and perform basic cleaning (section 11).

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## 6. Preparing the device

### **Attention - please note before assembling:**

- Check if any parts have been damaged in transit.
- For safe and trouble-free operation, only use manufacturer's original accessories and spares.
- Never use any aftermarket or damaged parts.

### 6.1 Charging the battery

With the charger unit connected to the mains, the instrument can be used immediately. For mains-independent use, the battery must be charged first.

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**Important:**

*To charge the battery, the device must be switched on and have a module attached.*

- Insert the plug of the charger unit into the charging socket.
- Connect the mains plug of the charger unit to the mains.
- After approx. 3 seconds, "Chg" will appear on the top right of the LC display, indicating the beginning of the charging process.

**Important user benefits:**

- The charging process automatically terminates when the battery is 100 % full (displayed in the dialog line). The battery cannot be overcharged.
- The instrument can be used normally during the charging process.

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
**Attention:**

Use only the original charger unit. The mains voltage must correspond to the voltage indicated on the recharger unit (maximum permissible deviation  $\pm 10\%$ ).

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## 6.2 How to Connect the Suction Tube

- Hold suction tube next to supply bottle with the retainer ring pointing towards the bottom of the bottle.

- Cut the tube at an angle at the height of the lower edge of the rim on the reservoir bottle.

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

- Attach suction tube securely.

- Next screw retainer ring onto suction valve.

### **6.3 How to Connect the Recirculation Tube**


- Firmly press recirculation tube in the designated recirculation outlet of valve block.


### **6.4 How to Connect the Discharge Tube**

- Discharge tube includes:

- tube guide
- retainer ring
- dispensing tube
- closure cap

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
- Firmly tighten retainer ring and verify correct positioning.

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
### **Please note:**

- Confirm correct positioning of the discharge tube. Only use original discharge tube (refer to accessories and replacement parts). Never install a deformed or defective discharge tube.

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### **6.5 Secure Device on Bottle**

- Screw device either with or without adapter onto reservoir bottle.

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### **6.6 Align Front of Device with Bottle Label**

- Rotate device to desired position of bottle label

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## **6.7 Align the discharge tube to point into the desired direction**


### **Warning!**

- Always follow all safety instructions and avoid prohibited and restricted applications.
- Always wear protective clothing, glasses and gloves.

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## **7. Start operation**

- Attach an electronic module (Dispenser or Titration) to the device prior to start operation.

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
### **Please note:**

- Always follow safety instructions (refer to section 1) as well as application prohibitions and restrictions (refer to sections 4 and 5).
- Only use the device in such a way that the discharge tube is never directed towards the user or other people. Therefore position the discharge tube with valve block to the right or left away from your body. Then align the operation module in a way that working with the keys is conveniently possible.

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

### **7.1 Remove Closure Cap**

- Place receiving container underneath discharge tube.
- Remove closure cap as shown in the figure.

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
**Please note:** Reagent may drip during this process!

- Slide closure cap back as shown in the figure.
- During operation the closure cap should rest in the safe storage position.

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
## 7.2 How to Prime the Device

Set the recirculation valve control to point toward the back of the device and press the PRIMING key in the even that air bubbles are occurring in the cylinder chamber e.g. after screwing the device onto a new reagent bottle. After pressing the key the piston moves approximately 25 mm

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
up, switches over and moves downward again. This process is repeated as long as the key is pressed. The system finishes this process after the key is released. The piston moves to the basic position in the event that it had stopped indefinitely (e.g. after pressing the STOP key). If air bubbles are no longer visible in the view glass, the recirculation valve control should be set to dispensing/titration (control points toward the discharge tube) and the actual dispensing/titrating may begin.

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### 8.1 Darstellung Elektronikmodul "dispenser"

## 8. Dispensing with the electronic module "Dispenser"

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
## 8.2 Function keys "Dispenser"

### I STOP

*Pressing this key stops the piston immediately. The LC display shows the residual volume that still has to be dispensed.*


*Pressing the DISPENSE key releases the discharge of this residual volume.*

### II CLEAR

*Pressing this key sets the lines displaying the dialog and the volume back to the starting point ("0").*

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
### III PRIMING

After pressing the key the piston moves approximately 25 mm up, switches over and moves downward again. The piston moves to the basic position in the event that it had stopped indefinitely (e.g. after pressing the STOP key). This process is repeated as long as the key is pressed. The system finishes this process after the key is released.

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
### IV REPEATED DISPENSING

The dialog-line in the LC display shows "t->000" after the REPEATED DISPENSING key is pressed. The desired number of repetitions is keyed in with the 3 operational keys on the left (100, 10 and 1

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## **V    PAUSE**


key). The entered value is either saved by pressing the **REPEATED DISPENSING** key again, or automatically after 7 seconds. Therefore up to 999 dispensing repetitions are possible. After repeated dispensing is programmed, the pause between each dispensing step may be set.

The dialog-line in the LC display shows "t->000" after the **PAUSE** key is pressed.

The desired pause is keyed in with the 3 operational keys on the left (100, 10 and 1 key). The entered value is either saved by pressing the **PAUSE** key again, or automatically after 7 seconds. Therefore pause time of up to 999 seconds.

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
## **VI DISPENSE**

This key starts either the single dispensing step or the set dispensing repetitions. The keyed-in volume is aspirated and dispensed immediately after this key is pressed.

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
### **8.3 Operational keys "Dispenser"**

The values

- I Volume
  - II Repetitions
  - III Pauses
- are set with these keys.

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## I VOLUME

Volumes between 1 to 500 ml (50 ml device) or 0,1 to 200 ml (20 ml device) are set by means of the respective operational key.

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
### **Please note:**

The LC display indicates the currently dispensed volume after the pistons has traveled down if the volume exceeds 50 ml (50 ml device) and respectively 20 ml (20 ml device). While the piston moves upward again, the keyed-in total volume and the volume dispensed so far are displayed alternately.

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
## II DISPENSING REPETITIONS


The number of desired dispensing repetitions is keyed-in with the 3 operational keys on the left in combination with the function key REPEATED DISPENSING.

Repetitions of up to 999 times are possible.

- Press the REPEATED DISPENSING key
- The dialog line display turns from "n=001" into "n->000"
- Enter the number of desired repetitions

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- The entered value is either saved by pressing the REPEATED DISPENSING key again, or automatically after 7 seconds.

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- The dialog line display turns from "n->000" into "n=xxx" (xxx = set value)

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### III PAUSE

The desired pause time between each dispensing step (only for repeated dispensing) is keyed-in with the 3 operational keys on the left in combination with the function key PAUSE. Pauses of up to 999 seconds are possible.

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- Press the PAUSE key
- The dialog line display turns from "t=000" into "t->000"
- Enter the desired pause

100	Set the 100-digit
10	Set the 10-digit
1	Set the 1-digit

- The entered value is either saved by pressing the PAUSE key again, or automatically after 7 seconds.

- The dialog line display turns from "t->000" into "t=xxx" (xxx = set value)

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
## 9.2 Function keys "Titration"

### I STOP

- STOP function during filling or priming:  
Pressing this key stops the piston immediately. The user is further guided with the dialog in the LC display.

- STOP function during fast titration with pre-selected volume:  
After pressing this key, the LC display shows the volume that has already been delivered.

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

## **III PRIMING**

*After pressing the key the piston moves approximately 25 mm up, switches over and moves downward again. The piston moves to the basic position in the event that it had stopped indefinitely (e.g. after pressing the STOP key). This process is repeated as long as the key is pressed. The system finishes this process after the key is released.*

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## **IV FILL**

### *- FILL without volume pre-selection*

*The glass cylinder is filled with the maximum volume 50 ml (50 ml device) respectively 20 ml (20 ml device) after pressing the key.*

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
### *- FILL with volume pre-selection*

*The volume (ranges between 2-50 ml), set with the operational keys (see section 6.3) is aspirated into the glass cylinder after pressing the key.*

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## **V DROPWISE TITRATION**

*Pressing the key briefly delivers volume dropwise. The dropwise titration is accelerated if*

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
pressed continuous. The delivered titration volume is permanently indicated in the LC display. The smallest possible titration volume is 0,01 ml.

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## VI FAST TITRATION

### 1. FAST TITRATION without volume pre-selection

Titration volume is delivered after pressing the key. The delivered titration volume is permanently indicated in the LC display.

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
### **Please note:**

In order not to titrate too much volume it is advisable to switch over to dropwise titration before the point of color change is achieved.

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## 2. FAST TITRATION with volume pre-selection

The volume, set with the operational keys is delivered after pressing the key.

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### **Please note:**

The glass cylinder is automatically refilled in the event that the amount of fluid in it is less than the set volume. The volume that has already been titrated is added up automatically.

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### **Please note:**

The final 10% of the pre-selected volume is automatically delivered dropwise. Once the delivery is completed, the titrated volume is indicated in the LC display.

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
### **Important advantage for the user:**

The ongoing titration is stopped abruptly after the titration point is achieved, if pressing the STOP key. The delivered titration volume is indicated in the LC display. Further proceeding with the keys FAST TITRATION or DROPWISE TITRATION is of course possible anytime.

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## **9.3 Operational keys "Titration"**

The values I Volume pre-selection during filling  
II Volume pre-selection for delivery with the function TITRATION WITH VOLUME PRE-SELECTION are set with these keys.

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
**Please note:**

Volumes between 1 to 50 ml (50 ml device) respectively 0,1 to 20 ml (20 ml device) can be set.

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
*I Fill variably*

The desired amount of fluid is keyed-in with the 4 operational keys 10, 1, 0,1 and 0,01 in combination with the functional key FILL. Pressing the FAST TITRATION key releases the delivery of the pre-selected volume.

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*II Volume pre-selection for delivery with the function TITRATION WITH VOLUME PRE-SELECTION*

*The desired amount of fluid is keyed-in with the 4 operational keys 10, 1, 0,1 and 0,01 in combination with the functional key FAST TITRATION. Pressing the FAST TITRATION key releases the delivery of the pre-selected volume.*

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
## **10. Special functions**

### **10.1 Guide language selection**

*The user can set the dialog line of the LC display by means of the operational keys to 4 different languages:*

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
- German
  - English
  - French
  - Spanish
- Selecting the guide language:*

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

- Switch off the device
- Press the operational key "0,01" and  
continue pressing it.*

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- Switch on the device and still continue  
pressing the operational key for approxi-  
mately 5 seconds.

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

- The guide language is set with the operational keys "10" German, "1" English, "0,1" French and "0,01" Spanish

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- The guide language set in this manner is saved.

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
The default setting of the device considers a reference temperature of 20°C with the cubical coefficient of expansion for borosilicate glass (DIN ISO 4787) and is calibrated "to deliver" (EX 20°C). Water of the "quality 3" in accord with DIN ISO 3696 is used for calibration. The user may exactly calibrate the device

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for any specific application (e.g. highly viscous fluids) himself, by means of the QUICK-CAL function. 50 ml (50 ml device) or respectively 20 ml (20 ml device) of the particular fluid are weighed on a sufficiently precise balance and the real volume is calculated under consideration

of the fluid's density. After that the device is switched off, the operational key "0,1" is pressed and continuously kept pressed while the device is switched on again. Still continue pressing the key for another 5 seconds.

QUICK-CAL appears in the display. The new adjustment is saved if the CLEAR key is pressed after the default setting 50.000  $\mu\text{L}$  (50 ml device) or respectively 20.000  $\mu\text{L}$  (20 ml device) is replaced with the calculated "real" volume.

Any volume is corrected automatically in accord with the performed calibration. A small "c" instead of the

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*decimal point refers to the differing calibration after any dispensing/titration.*

*Using the same procedure enables reset the default setting. In the event the value is reset to 50.000 µl (50 ml device) or respectively 20.000 µl (20 ml device), the decimal point appears again instead of the small "c".*

### **10.3 Adaptation to the basic unit**

*The dialog line of the operation module has to be adapted to the respectively used basic unit (20 ml or 50 ml) by means of the operational keys:*

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- Press the operational key "1" and continue pressing it.
- Switch on the device and still continue pressing the operational key for approximately 5 seconds.
- The module is programmed with the operational key 1 for the 20 ml device or respectively the operational key 10 for the 50 ml device. The current adjustment is indicated briefly any time the device is switched on.

#### **10.4 Adaptation of the dispensing speed**

The user can adapt the dispensing speed by means of the operational keys as follows:

- Switch off the device
- Press the operational key "10" and continue pressing it.

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- Switch on the device and still continue pressing the operational key for approximately 5 seconds.

Now the operational key 0,1 (slow), the operational key 1 (medium) or the operational key 10 (fast) are used to adapt the dispensing speed.

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
## **10.5 Adaptation of the suction speed**


The user can adapt the suction speed by means of the operational keys as follows:

- Switch off the device  
- Press the operational key "100" and continue pressing it.  
- Switch on the device and still continue

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pressing the operational key for approximately 5 seconds.

Now the operational key 0,1 (slow), the operational key 1 (medium) or the operational key 10 (fast) are used to adapt the suction speed.

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## 11. Cleaning

To warrant the lifetime of your device, it must be cleaned on a regular basis.

The following situations require the device to be cleaned:

- In the event the piston is jammed
- After titration of reagents whose application is restricted (refer to section 4+5)

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- Prior to change of reagent
- Prior to sterilization
- Prior to exchange of suction, discharge and recirculation valves
- Prior to longer storage periods

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
### **Please note:**

- Before the attempt is made to clean the device, put on protective clothing, glasses and gloves.
- Safety instructions (refer to section 1) must be observed at all times.
- Suction tube, discharge tube, cylinder and valves are filled with fluids. Direct all openings away from your body and other persons.

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## 11.1 Maintenance / Cleaning

- Complete dispensing/titration and put closure cap in place.
- Place the device, which is attached to a reservoir bottle into the cleaning vessel.
- Unscrew and remove device from reservoir bottle wearing suitable hand protection. (Suction tube should be out of the fluid in the bottle).

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

### **Please note:**


Fluid may spill out of the suction tube.

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
*Hold discharge tube above the opening of the reservoir bottle.*

- Remove closure cap and park in storage position.
- Switch on the device
- Press the **PRIME** key as long as no more fluid is discharged.

### **11.2 Basic cleaning by means of the priming-key**

- Dip suction tube into a suitable cleaning solution and rinse the device thoroughly by pressing the **PRIME** key. Keeping this key pressed will repeat this procedure

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**Caution:**

- Remaining fluids may splatter!
- Always direct the spout away from your body!

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
- Rinse the device with a suitable solvent (e.g. distilled water or acetone) and repeat the procedure by keeping the PRIME key pressed.

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- Remove suction tube from solvent and empty tube and device completely through pressing the PRIME key as long as no more fluid is discharged.

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
- Remove the following parts from the device:

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
- Remove recirculation tube

- Loosen suction tube retainer ring and remove suction tube. Place these parts in a suitable cleaning solution.

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### 11.3 Thorough Cleaning

- First perform cleaning as described in section 11.1 above.
- Then remove valve block, cylinder and piston:
- Rotate upper section of the device to the right until the movement stops (mark is lined up with the right edge of the chrome ring), then pull out upper section 6 mm and turn to the left until movement stops. Remove valve block with cylinder from upper section by pulling downward.

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#### **Please note:**

If necessary, piston and toothed rack can be separated by sliding the piston off the rack.

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
- Check piston for damage on the sealing lips (exchange if needed), then place piston into cleaning solution as described in section 11.1.


- Clean cylinder with brush and cleaning solution, then rinse with rinsing solution.

- Reassemble piston, toothed rack, valve unit and cylinder.

- Slide piston onto toothed rack.

- The piston may be moved back around 2 mm by pressing the DISPENSE or respectively the FAST TITRATION key provided that the device is switched on and a module is attached. This procedure may become advisable in the event that the glass cylinder and valve block will not go back into their place completely while reassembling the device.

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- Place upper section of the device onto the valve block with cylinder. The mark on the upper part must align with the left edge of the marking on the chrome ring. Insert upper section 6 mm into valve block with cylinder until stop, then rotate to the right until movement stops and push down another 6 mm. Move discharge spout to the desired working position.

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### **Please note:**

Carefully insert piston into the glass cylinder so that sealing lip will not be damaged. A damaged sealing lip will cause leakage problems!

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- Verify proper functioning of the device.

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## 12. Sterilization

All parts of the device which come in contact with fluid can be steam-sterilized up to 121 °C / 2 bar and sustain a sterilization time according to DIN te = 20 min.

### **Please note:**

Only the following components can be steam-sterilized:

- Valve block with cylinder
- Suction tube
- Discharge tube
- Recirculation tube
- Piston unit (very important: insert into glass cylinder prior to sterilization) Never steam-sterilize the upper section of the device or the electronic module! Refer to Individual Components section 2)

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## **12.1 Preparations for Sterilization**

- The device must be cleaned (refer to section 12) prior to sterilization.

## **12.2 Disassembly**

### **Please note:**

*The device has to be in the basic position (piston is all the way down) to disassemble it. The device has to be switched off for disassembly.*




- In addition to the parts listed under 11.1 all other parts of the device which come in contact with fluid must be disassembled and cleaned as described under 11.2.

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

### 12.3 How to Loosen Suction, Discharge and Return Valve

- Remove suction tube
- Remove discharge unit
- Loosen suction valve by turning it 180 degrees counter-clockwise with enclosed wrench.
- Loosen discharge valve by turning it 180 degrees counter-clockwise with enclosed wrench.
- Grasp cover at its sides and remove.

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- Loosen recirculation valve by turning it 180 degrees counter-clockwise with enclosed wrench.

- Loosen nut 180 degrees counter-clockwise and press on nut end of the stopcock to loosen stopcock.

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### *Remove the piston*

- Press threaded bar into the direction of the recirculation control and slide the piston into the other direction with the other hand.

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**Important:**

The piston has to be cleaned in a suitable solution prior to sterilization. The insert the piston into the glass cylinder, so that the coupling still rises from the rim of the glass cylinder.

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
**12.4 How to Perform a Sterilization**

Perform steam-sterilization of all parts, coming in contact with fluid (up to 121 °C/2 bar, sterilization time according to DIN  $t_e$ =20 min)

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
### **Please note:**

- Allow device to cool down to room temperature, then begin reassembly.
- Check all parts for damage and replace if needed.

## **12.5 Reassembly**

- Reassemble piston, valve block and cylinder:

- Place upper section of the device onto the valve block with cylinder. The mark on the upper part must align with the left edge of the marking on the chrome ring. Insert upper section 6 mm into valve block with cylinder until stop, then rotate to the right until movement stops and push down another 6 mm. Move dischar-

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ge spout to the desired working position.


- The piston may be moved back around 2 mm by pressing the DISPENSE or respectively the FAST TITRATION key provided that the device is switched on and a module is attached. This procedure may become advisable in the event that the glass cylinder and valve block will not go back into their place completely while reassembling the device.

**Please note:**

Insert the piston carefully into the glass cylinder. Otherwise its sealing lips may be damaged. This would cause leakage of the system.

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## 9.6 How to Secure Suction, Discharge and Recirculation Valve

- Tighten suction valve clockwise with wrench.
- Tighten discharge valve clockwise with wrench.
- Tighten recirculation valve clockwise with wrench.
- Secure nut on valve shaft clockwise by hand.
- Position cap.

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## **12.7 Functional test**


- Attach electronic module and switch on the device.
- Functional test with the PRIME key. The device must perform the priming procedure (also possible without fluid).
- In case of malfunction refer to "Troubleshooting" (section 17) or get in touch with the manufacturer.

## **13. How to Exchange Valves**

Leaking or damaged suction, discharge and recirculation valves must be exchanged immediately (refer to accessories and replacement parts in section 19).

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

### 13.1 How to Exchange the Suction Valve

- Clean device (refer to section 11).
- Disconnect suction tube.
- Turn device upside down. Otherwise valve ball may fall out during the exchange.
- Remove suction valve including valve ball with the help of the wrench.
- Insert ball and install new suction valve. Tighten with wrench.
- Reattach suction tube snugly.
- Next screw retainer ring onto suction valve (refer to section 6.2).

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### 13.2 How to Exchange Discharge Valve

- Clean device (refer to section 11).
- Disconnect discharge tube (refer to section 14.2).
- Unscrew discharge valve with wrench and remove.
- Screw in new discharge valve and tighten with wrench.
- Position discharge tube tightly and secure retainer ring (refer to section 6.4). Check for proper positioning!

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### 13.3 How to Exchange the Recirculation Valve

- Unscrew recirculation valve with wrench and remove.

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*- Screw in new recirculation valve and tighten with wrench.*

## **14. Switching Operation Titration/ Dispensing**

*By changing of the electronic modules and the discharge tubes the basic unit can be set from a dispenser to become a burette and reverse.*

### **14.1 How to change the modules**

*By merely changing the electronic module in the basic unit, the device can either be set for dispensing or titration. Make sure that the electronic modules are removed and inserted straight and properly.*

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## **14.2 How to exchange the discharge tube**

*Remove the attached discharge tube (refer to section 14.2) and attach another discharge tube, that is suitable for the desired purpose (refer to section 6.4)*

## **15. Repairs at the Manufacturer**

- *Prior to shipping a defective device to the manufacturer for repair it must be carefully cleaned (refer to section 11).*
- *Devices will only be repaired if they have been properly cleaned and decontaminated.*
- *Without proper proof of decontamination, devices will be considered contaminated. They may be serviced at higher costs or returned unserved.*
- *The malfunctioning of the device should be described in detail.*

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- Please indicate which reagents have been used. (Important!)
- If the device has been applied other than for its intended use or if the use of force is evident all warranty claims are null and void.

## 16. How to Test the Correct Volume

The parameters listed under section 18 for accuracy (R%) and coefficient of variation (CV%) of a device can be tested as follows:

### 1. Weighing process (refer to DIN 8655-6):

- Fill the burette completely with de-ionized water and discharge the nominal volume.
- Weigh discharge quantity on analytic scale.
- Convert quantity into volume.

**Note:** Temperature, air pressure and material must be taken into consideration.

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## **2. Repeated weighing process**

- The weighing process described above should be repeated at least five times to obtain reliable test results for the accuracy as well as coefficient of variation of a device.

## **3. Computation of accuracy and coefficient of variation**

- Accuracy (R%) and the coefficient of variation (CV%) of the device must be computed according to the formulas of statistical quality control.

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## 17. Trouble shooting

<i>Symptom</i>	<i>Possible cause</i>	<i>Solution</i>
<i>Fluid leaks at the gap between the valve block and the valves</i>	<i>- Valves loose or damaged</i>	<i>- Tighten valves with the wrench or replace if necessary (see point 13) Note: After sterilization (see point 12), suction/discharge /recirculation valves need to be tightened.</i>
<i>Fluid is leaking between the discharge tube and the valve</i>	<i>- Discharge tube loose or damaged - Wrong discharge tube used</i>	<i>- Mount the original discharge tube correctly (see point 6.4). - Use only the original discharge tube.</i>
<i>Device aspirates air</i>	<i>- Device not primed - Suction tube loose - Suction tube or retainer ring damaged</i>	<i>- Prime the device (see point 7.2). - Tighten suction tube (see point 6.2). - Replace suction tube (see accessories - spare parts point 19).</i>
<i>Fluid leaks at the gap</i>	<i>- Recirculation valve loose</i>	<i>- Tighten recirculation valve with the wrench</i>

<i>between the valve block and the recirculation valve</i>	<i>or damaged</i>	<i>or replace if necessary (see point 13.3)</i>
<i>Fluid leaks between stopcock shaft and valve block only perceptible when cover is removed</i>	<i>- Stopcock shaft loose</i>	<i>- Tighten nut on stopcock shaft</i>
<i>Delivered volume incorrect</i>	<i>- Discharge tube loose or damaged</i> <i>- Suction tube loose or damaged</i> <i>- Suction/discharge/recirculation-valves loose or damaged</i>	<i>- Mount original discharge tube correctly (see point 6.4).</i> <i>- Tighten suction tube (see point 6.2).</i> <i>- Tighten valves with the wrench.</i> <i>- Clean the device (see point 11).</i>
<i>Piston does not move up</i>	<i>- Jammed in bottom position</i>	<i>- Slightly turn valve block to the left or right, while holding the device's upper section</i>
<i>Fluid has passed the piston</i>	<i>- Damaged piston sealing lip</i>	<i>- Exchange piston cleaning/sterilization</i>
<i>Air bubbles visible in the glass cylinder</i>	<i>- Device not primed properly</i> <i>- Suction tube or retainer ring</i>	<i>- Prime the device (section 7.2)</i> <i>- Tighten retainer ring (section 6.2), loose damaged or deformed exchange original suction tube (see spares/accessories, section 19)</i>

<i>Volume incorrect</i>	<i>- Discharge tube loose or damaged</i>	<i>- Properly mount discharge tube (section 6.4)</i>
<i>Device does not react</i>	<i>- batteries are empty, CHG appears in the display</i>	<i>- Insert plug of the charger unit into the port of the device</i>
<i>Device does not charge</i>	<i>- Mains without current</i> <i>- Improper contact charger unit - mains socket</i> <i>- Electronic module not attached to device</i> <i>- Charger unit not plugged properly into device</i> <i>- damaged charger unit</i>	<i>- Check contact or use different mains socket</i>  <i>- Attach electronic module to device</i>  <i>- Plug charger unit properly into the port of the device</i> <i>- Send charger unit to manufacturer for test</i>
<i>No display after the device is switched on, despite fully charged batteries</i>	<i>- no contact between electronic module and device</i>	<i>- Remove electronic module, check pins and insert again. Assure proper contact by pressing onto the module's edges. Send electronic module to manufacturer for test</i>
<i>Piston moves stiffly</i>	<i>- No proper coupling of device</i>	<i>- refer to page 63, section 12.5 and threaded bar</i>
<i>Piston not in bottom position (not visible)</i>	<i>- incorrect operation actuated the on/off switch</i>	<i>- Switch the device off/on. Press the prime key. Warning: Fluid will be discharged.</i>

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**General advice:**

- Briefly switch the device off and on again, in the event of incorrect dialog display or other malfunctions of the device. This procedure checks the system software and the device's components (RESET function).
- The battery capacity (in %) in the LC display is only an approximate indication of the real battery capacity.

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**18. Technical data****Basic unit:**

Type:

Modular, micro-processor controlled piston-system with reagent recirculation for mains-independent battery operation.

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*1 ml - 500 ml (maximum piston displacement 50 ml, adding automatically). Resolution: 7 µl*

*Volume accuracy:*

*Dispensing A\* 0,6%, CV\* 0,1%*

*Titration A\* 0,2%, CV\* 0,1%*

*\* Stated values for accuracy (A) and coefficient of variation (CV) refer to the maximum volume 20 ml or respectively 50 ml at an equal temperature (20°C) of device, ambient and fluid (de-ionized water).*

*Dispensing rate: Adjustable in 3 stages*

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#### *QUICK-CAL:*

*Calibration for other temperatures and fluids possible*

#### *Power supply:*

*Built-in high performance battery pack or with mains adapter 15V 850 mA*

#### *Battery operation:*

*Capacity per charge approximately 15.000 ml. Battery capacity indicated in the LC display, charging time around 2,5 h, intelligent technology rules overcharging out.*

#### *Material:*

*Case made of PA6, all parts getting in touch with fluid made of chemically resistant material (ECTFE, DURAN, Hastelloy, FEP, PTFE)*

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*GS-geprüft:*

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*Cleaning/Sterilization:*

*Motorized rinsing function by pressing a key, decoupling of dry and wet parts without tools, all parts getting in touch with reagent are steam-sterilizable.*

*Priming:*

*Automatically by pressing a key*

*Dimensions:*

*350 mm x 170 mm x 95 mm / 900 g (with electronic module)*

*Certification:*

*Come up to meet ISO 9000 certification with its individual test protocol with serial number.*

*GS-approved:*

*The device is GS-approved by the "TÜV-Süd"*

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*Conformity:*

*Conformity certified in accord with the German "Eichordnung" (Appendix 12)*

*CE:*

*This sign confirms that the device was manufactured and tested in accord with the existing EU-guidelines.*

*Warranty:*

*12 months warranty*

### **Electronic modules:**

*Module type:*

*Specific electronic modules for dispensing or respectively titration  
The electronic module includes the complete control electronic (micro-processor, storage chip etc.)*

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*Key pad:*

*Foil key pad with integrated 2-line LC display for entering values and commands.*

*User guidance:*

*Clear functional reference of any key, dialog guidance, guide language adjustable (D, GB, F, E)*

**Electronic module "Dispenser":**

*Single dispensing:*

*Volume from 1 ml to 500 ml (50 ml device)  
or respectively 0,1 ml to 200 ml (20 ml device)*

*Time-controlled serial dispensing:*

*Volume from 1 ml to 500 ml (50 ml device)  
or respectively 0,1 ml to 200 ml (20 ml device)  
Number of repetitions (1 to 999), Pauses (1 s to 999 s)*

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### **Electronic module "Titration":**

*Volume pre-selection:*

*Up to 50 ml (50 ml device) or respectively 20 ml (20 ml device)*

*Fast titration:*

*Speed adjustable in 3 stages (see special functions)*

*Dropwise titration: key with 1 drop discharge per actuation, accelerating when key is pressed continuously*

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### *Discharge unit*

#### *Code-Nr.*

931 30 25	Dispenser discharge, completely guided tube 115 mm
931 30 15	Titration discharge, completely guided tube, pointed 115 mm
956 30 04	Titration, variabel, with tube holder, approx. 400 mm
931 30 27	spiral, working length approx 1 m

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*Discharge valve*

*Code-Nr.*

931 60 05	valve seat/ball from glass, spring from Hastelloy
931 60 29	valve seat/ball from glass, spring from platinum/ iridium

*Suction valve*

*Code-Nr.*

931 65 05	valve seat/ball from glass, spring from Hastelloy
931 65 29	valve seat/ball from glass, spring from platinum/ iridium

*Recirculation valve*

*Code-Nr.*

931 66 25	valve seat/ball from glass, spring from Hastelloy
931 66 29	valve seat/ball from glass, spring from platinum/ iridium

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<i>Code-Nr.</i>	<i>Volt</i>	
956 75 01	230 V	EURO
956 76 01	240 V	UK
956 85 01	120 V	US

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