**Eppendorf Varipette® 4720**  
**Bedienungsanleitung**  
**Instruction Manual**  
**Mode d'emploi**  
**Istruzioni d'impiego**  
**Manual de Instrucciones**

**eppendorf**

Eppendorf AG  
22331 Hamburg · Germany  
Phone: +49 40 538 01 0  
Fax: +49 40 538 01 556  
E-mail: [eppendorf@eppendorf.com](mailto:eppendorf@eppendorf.com)  
Internet: [www.eppendorf.com](http://www.eppendorf.com)  
Application Hotline:  
Phone: +49 180 366 67 89  
E-mail:  
[application-hotline@eppendorf.com](mailto:application-hotline@eppendorf.com)

B 4720 900.012-05/0104  
Printed in Germany

**eppendorf**  
eppendorf® is a registered trademark

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13

<p><b>Eppendorf Varipette® 4720</b></p> <p>This instrument is part of our extensive liquid handling program.</p> <p>The Varipette 4720 has been manufactured and tested under strict quality controls to ensure problem-free operation. Please take a few minutes and read this instruction manual in order to familiarize yourself with the features of the Varipette.</p>	<p><b>1 Introduction</b></p> <p>The Eppendorf Varipette 4720 is a versatile, variable milliliter pipette. Volumes from 1.00 to 10.00 ml can be set in 10 µl increments by turning the volume selection dial.</p> <p>The Varipette can be used to pipette aqueous and organic solutions as well as liquids with a high viscosity or solutions whose density and vapor pressure vary considerably from those of water.</p> <p>The 4720 becomes a functional unit when a Varitip is attached. The Varitip consists of a piston, barrel and tip and holds the liquid to be dispensed.</p> <ul style="list-style-type: none"> <li>- The Varitip P (1.00 – 10.00 ml) with a short tip works according to the positive-displacement principle.</li> </ul>	<ul style="list-style-type: none"> <li>- The Varitip S with Maxitip is for pipetting in the range from 2.50 ml to 10.00 ml. The long Maxitip is designed for aspirating liquids from tall and long-necked vessels (diameter of tip = 10 mm, fits into flasks with standard ground joint 12/21). During pipetting, solution only enters the Maxitip, therefore, if the Maxitip is exchanged, samples and reagents can be dispensed in succession without carryover.</li> <li>- If liquids with a high vapor pressure or density are to be dispensed, the required volume setting can easily be determined using the graduated Maxitip G.</li> </ul>	<ul style="list-style-type: none"> <li>- The valve which can be inserted into the Maxitip enables more problem-free pipetting of aggressive or toxic solutions or liquids with radioactive or pathogenic components.</li> <li>- The Maxitip is made of polypropylene and autoclavable. The barrels of the Varitips are made of polypropylene and the pistons of polyethylene.</li> </ul> <p>The Pipette Stand 3115 is designed to hold the Varipette 4720 with the aid of an adapter.</p>	
---	--	--	---	---

14

	<p><b>1.1 Design principle</b></p> <p><b>The Varipette</b> (Fig. 2):</p> <p>(1) Volume selection dial</p> <p>(2) Pipetting lever Liquid is aspirated by pressing this lever down to the stop. It is dispensed by pressing the lever again.</p> <p>(3) Volume readout</p> <p>(4) Locking clamp Serves to clamp the piston firmly into the Varipette.</p>	<p><b>The Varitips</b> (Fig. 2):</p> <p>(5) Varitip P, 1.00 ml – 10.00 ml (= Positive displacement)</p> <p>(6) Varitip S, dispensing part, 2.50 ml – 10.00 ml (= Varitip System)</p> <p>(7) Varitip S, Maxitip. The Maxitip is attached to the dispensing part and is interchangeable.</p> <p>(8) Graduated Maxitip G. Serves to set the Varipette volumetrically.</p>		<p><b>Valve in Maxitip</b> (Fig. 3):</p> <p>Can be inserted either way round into the Maxitip. The valve closes the Maxitip and prevents dripping.</p>
--	---	--	--	--

15

<p><b>1.2. Range of application</b></p> <p>Varitips P and Varitips S complement each other. Their recommended applicational ranges are shown in the following table:</p> <table border="1"> <thead> <tr> <th></th> <th><b>Varitip P</b></th> <th><b>Varitips S</b></th> </tr> </thead> <tbody> <tr> <td>Aqueous solutions</td> <td>+</td> <td>+</td> </tr> <tr> <td>Liquid aspiration from tall vessels</td> <td></td> <td>+</td> </tr> <tr> <td>Organic liquids (alcohol, organic acids, etc.)</td> <td>+</td> <td>+</td> </tr> <tr> <td>Solvents with high vapor pressure, e. g. diethyl ether</td> <td></td> <td>+*</td> </tr> <tr> <td>Liquids with high density</td> <td>+</td> <td>+*</td> </tr> <tr> <td>Viscous liquids</td> <td>+</td> <td></td> </tr> </tbody> </table> <p>* after using Maxitip G</p>		<b>Varitip P</b>	<b>Varitips S</b>	Aqueous solutions	+	+	Liquid aspiration from tall vessels		+	Organic liquids (alcohol, organic acids, etc.)	+	+	Solvents with high vapor pressure, e. g. diethyl ether		+*	Liquids with high density	+	+*	Viscous liquids	+			<p><b>2 Operation</b></p> <p><b>2.1 Inserting and removing the Varitip</b></p> <p><b>Important:</b> When inserting and removing the Varitip, always hold down the pipetting lever (see. Fig. 4). During insertion, the piston of the Varitip must be pushed fully into the barrel.</p> <ul style="list-style-type: none"> <li>– Press pipetting lever downwards and hold (locking clamp in down position). See Fig. 4.</li> <li>– Raise locking clamp.</li> <li>– Insert Varitip P or dispensing part of the Varitip S until it clicks into position. A slight twisting motion may aid this procedure. To remove Varitip, pull out of the mounting.</li> </ul>	<p>During insertion of the Varitip, the slipping clutch of the piston slides between the grooves of the locking clamp. If necessary, push in properly.</p> <ul style="list-style-type: none"> <li>– Lower locking clamp and press closed.</li> <li>– Let pipetting lever slide back slowly.</li> <li>– For Varitip S: If desired, insert the valve into the Maxitip (either way round). Attach the Maxitip to the dispensing part.</li> </ul> <p>Always let the pipetting lever slide back slowly to the stop. Do not let it snap back.</p> <p><b>eppendorf</b></p>
	<b>Varitip P</b>	<b>Varitips S</b>																						
Aqueous solutions	+	+																						
Liquid aspiration from tall vessels		+																						
Organic liquids (alcohol, organic acids, etc.)	+	+																						
Solvents with high vapor pressure, e. g. diethyl ether		+*																						
Liquids with high density	+	+*																						
Viscous liquids	+																							

16

	<p><b>2.2 Volume setting</b> (Fig. 5)</p> <p>The volume to be pipetted is set by turning the volume selection dial. The digits of the volume readout must appear completely in the window.</p> <p>The four-digit volume readout shows the dispensing volume in ml with two places after the decimal point.</p> <p>Example: Desired volume = 8.75 ml Readout: 08.75</p>	<p><b>2.3 Pre-wetting</b></p> <p>For accurate pipetting, the Varitip P and the Maxitip must be wetted first with the liquid to be pipetted. To achieve this, aspirate and dispense liquid three times. The solution can then be pipetted accurately.</p>	<p><b>2.4 Pipetting with Varitip P</b></p> <p>During filling, the piston must be in the down position.</p> <p><b>Filling</b></p> <ul style="list-style-type: none"> <li>- Immerse the Varitip P into the liquid.</li> <li>- Slowly press down the pipetting lever completely to the stop.</li> <li>- Hold for 2 sec. Slowly release the pipetting lever all the way up.</li> </ul>	<p><b>Dispensing:</b></p> <ul style="list-style-type: none"> <li>- Liquid can be dispensed from the Varitip P freely or against the inside of the vessel.</li> <li>- Slowly press down the pipetting lever completely to the stop.</li> <li>- Any remaining droplets on the Varitip P should be left along the inside of the vessel.</li> <li>- Now let the pipetting lever slide back slowly.</li> </ul>
--	--	--	--	---

17

<p><b>Note:</b></p> <p>After filling, a small air bubble remains in the barrel of the Varitip P. It has <b>no</b> measurable influence on the accuracy and precision of dispensing.</p> <p>Any liquid remaining in the Varitip P does not belong to the measured volume.</p>			<p><b>2.5 Pipetting with Varitip S</b> (Fig. 6 + 7)</p> <p>The Varitip S is filled and emptied in the same way as the Varitip P. However, please observe the following additional points:</p> <ul style="list-style-type: none"> <li>- When filling, always hold the pipette vertically.</li> <li>- When dispensing, hold the Maxitip against the inside of the vessel.</li> </ul>	<p>If, after dispensing strong wetting liquids, liquid remains in the Maxitip, it can be dispensed by an additional stroke:</p> <ul style="list-style-type: none"> <li>- Let pipetting lever slide back.</li> <li>- Move the piston upwards by slightly pressing the pipetting lever.</li> <li>- Dispense liquid by pressing the pipetting lever again.</li> </ul> <p><b>eppendorf</b></p>
--	--	--	--	--

18

**Pipetting with Maxitip and valve**

The valve closes the Maxitip with its own weight (density = 1.78).

- Before filling, hold the Maxitip against the inside of the vessel so that the valve is pushed back once.
- During filling, the flowing liquid lifts the valve. When dispensing or moving the piston downwards, hold the Maxitip at an angle against the inside of the vessel so that the valve is pushed up (see Fig. 7).
- Hold the Varitip with filled Maxitip and valve vertically.

**2.6 Pipetting liquids with special physical characteristics**

Liquids whose density and vapor pressure differ considerably from those of water can be dispensed with the Varitip S. For this purpose, determine the required volume setting with the graduated Maxitip G (see Fig. 8).

- Attach Maxitip G. Set desired volume by turning the volume selection dial.
- Aspirate liquid. Hold down the pipetting lever.
- By adjusting the volume selection dial, lower or raise the liquid meniscus to the appropriate graduation of the Maxitip G.

A counterclockwise rotation raises the liquid. A clockwise rotation lowers the liquid.

- Read from the volume readout the value which should be set to obtain the desired volume.

The liquid can then be dispensed using a Maxitip with or without a valve.

**Note:**

When pipetting liquids with a high vapor pressure, take up and eject liquid until no more gas bubbles are pressed out of the tip.

It is advisable to use the Varitip P when dispensing viscous solutions.

**3 Safety precautions and operational limitations**

Proper operation of the Varipette constitutes no danger in itself. However, volume errors may occur if there is a difference between the temperature of the Varitip system and the liquid to be dispensed.

**4 Gravimetric control of accuracy and precision**

The basis for gravimetric testing is DIN regulation 12 650. Bidistilled water which has taken on the ambient temperature of 20 to 25 °C (reference temperature, see Sec. 7) must be pipetted. To calculate the volume from the weight of the pipetted liquid, divide its weight by the density of the water.

5 measurements are sufficient to check the accuracy and 15 to determine the precision.

Due to its higher inaccuracy, the Maxitip G should not be used for testing the Varipette (see Sec. 7).

If liquids with special physical characteristics are to be pipetted, the Varipette can be set gravimetrically instead of volumetrically (see Sec. 2.6) by weighing the liquid pipetted when the nominal value is set.

New set value = 2 x nominal value –

$$\frac{\text{Weight of pipetted liquid}}{\text{Density of pipetted liquid}}$$

Check this value gravimetrically and correct if necessary.

**eppendorf**

## 5 Troubleshooting

Error	Cause	Elimination
Impossible to press pipetting lever down fully.	Locking clamp not pressed fully down.	– Slide piston properly into piston, close locking clamp completely.
Piston cannot be moved with the pipetting lever.	Pipetting lever is not fully in down position.	– Press pipetting lever down to the stop.
Incorrect pipetting volume.	Pipetting lever was not fully in down position during liquid aspiration.	– Always press pipetting lever fully to the down position. This position varies according to the setting for the volume to be pipetted.
	Locking clamp not completely down.	– Insert piston properly. Press locking clamp completely down.

Error	Cause	Elimination
Incorrect pipetting volume.	Liquid with deviating density or high vapor pressure.	– Perform volumetric or gravimetric control (see Sec. 2.6 and 4).
Dripping Varitip.	Leaking piston.	– Replace Varitip.
	Maxitip is not attached firmly to dispensing part.	– Attach tip firmly.
Droplets in the Maxitip.	Liquid with high vapor pressure or low viscosity.	– Insert valve into Maxitip.
	Too much wetting of the plastic.	– Replace Maxitip.

21

## 6 Maintenance and service

### 6.1 Maintenance

Varipettes should be cleaned with warm water or isopropyl alcohol. Apart from this, the Varipette is completely maintenance-free.

Individual sections of the Varipette should never be unscrewed!

### 6.2 Service

If a problem cannot be solved with the aid of the suggestions mentioned here, return the Varipette to your Eppendorf distributor for repair or replacement.

## 7 Technical data

### Varipette 4720 + Varitips S

Nominal volume	Systematic error (Inaccuracy)	Random error (Imprecision; CV)
2.5 ml	± 1.0 %	≤ 0.2 %
5 ml	± 0.4 %	≤ 0.2 %
10 ml	± 0.3 %	≤ 0.2 %

### Varipette 4720 + Varitip P

Nominal volume	Systematic error (Inaccuracy)	Random error (Imprecision; CV)
1.0 ml	± 0.6 %	≤ 0.2 %
5 ml	± 0.5 %	≤ 0.1 %
10 ml	± 0.3 %	≤ 0.1 %

For other volume ranges, the limits of the next smaller volume are valid.

Inaccuracy of Varipette 4720 and Varitip S when volume set with Maxitip G: =  
 ± 2.0 % at 2.5 ml  
 ± 1.5 % at 5.0 and 10.0 ml.  
 For imprecision values, see above.

Liquid: dist. water.  
 Reference temperature: 20 – 25 °C, constant to ± 0.5 °C  
 Number of determinations: 10, in accordance with ISO 8655 using Original Eppendorf Varitips P and Varitips S without valve, liquid dispensed against the inside of the weighing vessel.

Technical specifications subject to change!

**eppendorf**

22

<p><b>8 Ordering Information</b></p> <p>Eppendorf Varipette® 4720          Delivery package includes:          1 Varitip P,          1 Varitip S (1 dispensing part + 1 Maxitip),          10 valves,          1 graduated Maxitip G</p> <p><b>Disposables and accessories</b></p> <p>Varitip P (set of 100)          Varitip S, starter kit          (100 Maxitips, 10 dispensing parts,          10 valves, 2 graduated Maxitips G)          Dispensing parts for Varitip S (set of 30)          Maxitips autoclavable, for Varitip S (set of 200)          Maxitip G, graduated (set of 5)</p>		<p>4720 000.011</p> <p>0030 048.130</p> <p>0030 050.525</p> <p>0030 050.533</p> <p>0030 050.541</p> <p>0030 050.576</p>	<p>Pipette carousel,          (incl. 6 pipette supports) 3115 000.003</p> <p>Holder (for 4720, 4780 and 4980/81) 4981 201.006</p> <p>Holder (with adhesive surface          for attachment to lab benches and wall) 4981 203.009</p> <p>Please only use the accessories recommended by Eppendorf. Using spare parts or disposables which we have not recommended can reduce the precision, accuracy and life of our pipettes. We do not honor any warranty or accept any responsibility for damage resulting from such action.</p>	
---	--	---	--	--

23

<p><b>9 Short Instructions</b></p> <p>These short instructions may only be used if the operator is thoroughly acquainted with the detailed instructions, otherwise errors may occur.</p>	<p><b>Inserting the Varitip:</b></p> <ul style="list-style-type: none"> <li>- Select Varitip P or Varitip S.</li> <li>- Insert piston fully into the Varitip.</li> <li>- Press pipetting lever downwards and hold (locking clamp in down position).</li> <li>- Insert Varitip.</li> <li>- Lower locking clamp.</li> <li>- Let pipetting lever slide back slowly.</li> <li>- Varitip S: Insert valve into the tip, if required. Attach Maxitip to the dispensing part. Hold tip against the inside of the vessel so that the valve is pushed back once.</li> </ul>	<p><b>Dispensing:</b></p> <ul style="list-style-type: none"> <li>- Set pipetting volume.</li> <li>- Immerse Varitip P/Maxitip into the liquid.</li> <li>- To aspirate liquid, press pipetting lever slowly down to the stop.</li> <li>- Let pipetting lever slide back.</li> <li>- To pipette liquid, slowly press pipetting lever down to the stop. During dispensing, hold Maxitip against the inside of the vessel. If a valve is used, hold Maxitip at an angle so that the valve is pushed up.</li> <li>- Let pipetting lever slide back slowly.</li> </ul>		
--	---	--	--	---

24