

NEW!



Transferpette® *pro*

A New Twist on Comfort and Control
Simple. Easy. Efficient

BRANDTECH.com

Transferpette® *pro* micropipette

The addition of the Transferpette® *pro* to the trusted Transferpette® family offers scientists even more choice in selecting the ideal pipette for their unique needs. Designed for repeatability and comfort, the Transferpette® *pro* is built to meet the demands of your daily lab routine.

The Transferpette® *pro* features a slim, ergonomic handle paired with a rotating finger rest— a unique design that adjusts to fit any hand position, left- or right-handed. Combined with two intuitive calibration systems— Easy Calibration for tool-free, permanent adjustment and User Adjustment for quick, reversible fine-tuning—the Transferpette® *pro* ensures precise performance, even with changing liquid types, temperatures, or style of tips.

Your benefits:

- + **Simple:** Comfortable pipetting with a slim profile and adjustable finger rest
- + **Easy:** Adapt quickly to changing conditions with new user adjustment
- + **Efficient:** Reliable, durable, and convenient

A New Twist on Comfort and Control



Every application has different requirements. The Transferpette® *pro* adapts perfectly to your individual needs with additional innovative features:

Adjustable finger rest



Adjustable to 9 positions for an ergonomic hand position

User Adjustment



Easy and temporary adjustment for special conditions such as temperature, tip shape, or liquid type

Counter protection



Protects the pipette by preventing over-rotation, ensuring lasting performance and reliability

Extended volume range



Now available in 300 µl, 2.5 ml, and 1250 µl nominal volume

Counter protection

Prevents damage from over-rotating beyond volume range

Rotatable finger rest

Simple adjustment without tools

User Adjustment

Temporary adjustment to changing conditions

Extended volume range

13 single-channel pipettes from 0.1 µl to 10 ml
6 multi-channel pipettes from 0.5 µl to 1250 µl



All product information can be found at [BRANDTECH.com/pro](https://www.brandtech.com/pro)

Transferpette® *pro* and Transferpette® *S* in comparison

	Transferpette® <i>pro</i>	Transferpette® <i>S</i>
Volume range single-channel (adjustable)	12 sizes: from 0.1 µl to 10 ml	10 sizes: from 0.1 µl to 10 ml
Volume range multi-channel	6 sizes: from 0.5 µl to 1250 µl	5 sizes: from 0.5 µl to 300 µl
Size of handle	Slim	Standard
Adjustable finger rest	✓	-
One-handed operation	✓	✓
Uses universal tips	✓	✓
Volume lock	✓	✓
Counter protection	✓	-
Easy Calibration Technology	✓	✓
Temporary User Adjustment	✓	-
Completely autoclavable (121 °C)	✓	✓
UV-resistant	✓	✓
Warranty for single-channel pipettes	3 years	3 years
Warranty for multi-channel pipettes	1 year	1 year

