

JANUARY 19 – JULY 30, 2026

Labnet LabValues 2026

Save big with special offers for laboratory equipment and consumables



LV-1.26

Labnet™ MultiGene™ OptiMax

The MultiGene OptiMax Thermal Cycler delivers advanced speed and features while providing premium performance at an affordable price.

The OptiMax comes with standard built in 2-step, temperature optimization, touchdown and time increment protocols that are easy to adjust to meet your cycling needs.

- Unique flexible programming with FlexTemp™ technology
- Fast ramp rates 5°C heating / 3.5°C cooling

Cat. No.	Description	List Price	Special Price
TC9610	MultiGene OptiMax with 96-well block, 120V	\$4,118.00	\$2,382.00

Labnet AccuSeal™ Semi-automated Plate Sealer

Offering complete versatility, the AccuSeal will accept a full range of plates for PCR, assay or storage applications, and can be used for standard height to deep well plates.

- Adjustable sealing temperature: 100°C to 190°C
- Precise temperature, timing, and pressure for consistent sealing

Cat. No.	Description	List Price	Special Price
PS1000	AccuSeal semi-automated plate sealer, includes adapters for standard and deep well microplates, 120V	\$4,800.00	\$2,300.00



LV-2.26



LV-3.26

Labnet Spectrafuge™ 16M Microcentrifuge

The Spectrafuge 16M is the perfect high speed compact microcentrifuge. Producing a maximum of 16,000 x g, it is the ideal unit for DNA applications. A brushless, maintenance-free motor provides extremely fast acceleration and deceleration of the 18-place rotor and makes the Spectrafuge 16M one of the quietest microcentrifuge in its class.

Cat. No.	Description	List Price	Special Price
C0160	Spectrafuge 16M microcentrifuge, standard gray, 120V	\$1,499.79	\$899.87
C0160-R	Spectrafuge 16M microcentrifuge, red, 120V	\$1,499.79	\$899.87

Labnet™ Spectrafuge™ 6C Centrifuge

The Spectrafuge 6C has been designed for quick production of PPP (platelet-poor plasma) or PRP (platelet-rich plasma) as well as other applications in the clinical and research laboratory.

The fixed angle rotor of the Spectrafuge 6C holds up to 6 x 15 mL conical tubes of maximum tube height of 120 mm, and with adapters will accommodate most common blood draw tubes. Variable speed allows for a variety of applications to be performed.

- Large digital displays for time and speed
- Small footprint
- Cold room compatible

Cat. No.	Description	List Price	Special Price
C0060	Spectrafuge 6C centrifuge, with 6 x 10/15 mL rotor, 120V	\$1,214.28	\$850.00



Labnet Mini Microcentrifuge

Simple and convenient to operate, the Labnet Mini Microcentrifuges feature greater capacity (8-place microcentrifuge tubes and 4-place PCR strip rotors), a quick-release rotor system, and electronic brake. Once samples are loaded and the lid is closed, the rotor rapidly accelerates to 6,000 rpm, making it an ideal centrifuge for bringing small droplets to the bottom of tubes for micro-filtrations or basic separations.

- Quick release rotor system – no tools needed to exchange rotors
- Electronic braking – allows for more rapid deceleration, less waiting

Cat. No.	Description	List Price	Special Price
C1601-G	Mini centrifuge (grey lid), 100-240V	\$230.41	\$161.29
C1601-P	Mini centrifuge (purple lid), 100-240V	\$230.41	\$161.29
C1601-R	Mini centrifuge (red lid), 100-240V	\$230.41	\$161.29

Labnet FastPette™ Pro Pipet Controller

The FastPette Pro is a lightweight and ergonomically designed automatic pipet controller ideal for liquid handling with 0.5 to 100 mL glass or plastic serological pipets. Conveniently positioned switches allow toggling of operation modes and aspiration speeds for high precision handling to suit different pipet volumes and liquid viscosities. Additionally, aspirating and dispensing speed is controlled by two comfortable finger triggers. Liquids can be dispensed by gravity or by blow out using the motorized pump.

- Triggers convenient for one or two finger control of speed and mode
- LCD display for convenient visualization of speed, modes, and battery status
- Li-ion battery for fast charging and longer life

Cat. No.	Description	List Price	Special Price
P2004	FastPette Pro is supplied with a universal power supply*, an additional autoclavable 0.2 µm filter, and a wall hanger	\$449.12	\$179.65

*Universal power supply comes with US, EU, UK and AU plugs.



Labnet™ FastPette™ V2 Pipet Controller

The FastPette V2 is a battery powered pipet controller ideal for liquid handling with 0.5 to 100 mL glass or plastic serological pipets. Designed with the user in mind, the FastPette has a lightweight handle with smooth finger triggers and ergonomically located switches.

- Works with pipets from 0.5 mL to 100 mL
- Cordless work for up to 8 continuous hours
- Autoclavable pipet holder

Cat. No.	Description	List Price	Special Price
P2000	FastPette V2 Pipet Controller is supplied with a universal power supply*, wall hanger, bench stand, and 0.2 µm PTFE filter	\$358.06	\$143.22

*Universal power supply comes with US, EU, UK and AU plugs.



Laboratory Plastics and Consumables

Cat. No.	Description	List Price	Special Price
K213	Drying Rack	\$367.07	\$183.54
P9803	Krystal plate, 96 x 350 µL, black, tissue culture-treated, individually wrapped with lid, 100 pcs	\$1,107.44	\$775.21
R8300-BLK	100-place microtube box with hinged lid, opaque black for light sensitive samples	\$13.39	\$6.70

LV-8.26

Labnet Excel™ Electronic Pipettors

The Excel is an ergonomic, intuitive electronic pipettor for use in most laboratory applications. Its high precision stepping motor controls piston action, resulting in exceptional accuracy and precision.

- Factory calibrated to strict ISO8655 standard
- Large LCD screen and simple interface
- Ergonomically designed handle and finger rest

Cat. No.	Description	List Price	Special Price
P3600L-10	Excel. single-channel electronic pipettor. 0.5-10 µL	\$406.20	\$243.72
P3600L-20	Excel. single-channel electronic pipettor. 2-20 µL	\$406.20	\$243.72
P3600L-200	Excel. single-channel electronic pipettor. 20-200 µL	\$406.20	\$243.72
P3600L-1200	Excel. single-channel electronic pipettor. 100-1200 µL	\$406.20	\$243.72
P3608L-10	Excel. 8-channel electronic pipettor. 0.5-10 µL	\$829.44	\$497.66
P3608L-20	Excel. 8-channel electronic pipettor. 2-20 µL	\$829.44	\$497.66
P3608L-200	Excel. 8-channel electronic pipettor. 20-200 µL	\$829.44	\$497.66
P3608L-1200	Excel. 8-channel electronic pipettor. 100-1200 µL	\$940.27	\$564.16
P3612L-10	Excel. 12-channel electronic pipettor. 0.5-10 µL	\$968.18	\$580.91
P3612L-20	Excel. 12-channel electronic pipettor. 2-20 µL	\$968.18	\$580.91
P3612L-200	Excel. 12-channel electronic pipettor. 20-200 µL	\$968.18	\$580.91
P3612L-1200	Excel. 12-channel electronic pipettor. 100-1200 µL	\$1,106.31	\$663.79



LV-9.26



WHILE SUPPLY LASTS

Labnet™ Labpette™ FX Fixed Volume Pipettors

Labnet Labpette FX Fixed Volume Pipettors are ideal for applications where specific volumes are used repeatedly and sterilization by autoclaving is required.

The pipettors feature an ergonomically designed body and finger rest for comfortable grip even during extended pipetting sessions. The smooth plunger movement and unique twist-eject tip removal system require minimal force and help to reduce the risk of repetitive strain injuries.

- ▶ Factory calibrated to strict ISO 8655 standards
- ▶ Seventeen fixed volumes from 2 to 1,000 µL
- ▶ Ergonomically designed handle and finger rest

Cat. No.	Description	List Price	Special Price
P3900-20	Labpette FX Fixed Volume Pipettor, 20 µL	\$117.55	\$47.02
P3900-50	Labpette FX Fixed Volume Pipettor, 50 µL	\$118.72	\$47.49
P3900-250	Labpette FX Fixed Volume Pipettor, 250 µL	\$118.72	\$47.49
P3900-300	Labpette FX Fixed Volume Pipettor, 300 µL	\$118.72	\$47.49



LV-10.26

WHILE SUPPLY LASTS

Labnet BioPette™ A Pipettors

Available in single-, eight- and twelve-channel configurations, all BioPette A pipettors feature smooth plunger movement with low pipetting forces.

- ▶ Factory calibrated to strict ISO 8655 standards
- ▶ Fully autoclavable and UV resistant
- ▶ Designed for universal fit with all common brands of pipet tips



LV-11.26

Cat. No.	Description	List Price	Special Price
P3960-10A	BioPette A, single-channel pipettor, 0.5-10 µL	\$241.01	\$144.61
P3960-20A	BioPette A, single-channel pipettor, 2-20 µL	\$241.01	\$144.61
P3960-100A	BioPette A, single-channel pipettor, 10-100 µL	\$241.01	\$144.61
P3960-200A	BioPette A, single-channel pipettor, 20-200 µL	\$241.01	\$144.61
P3960-1000A	BioPette A, single-channel pipettor, 100-1,000 µL	\$241.01	\$144.61
P3960-5000A	BioPette A, single-channel pipettor, 1,000-5,000 µL	\$277.36	\$166.42
P3960-10000A	BioPette A, single-channel pipettor, 1,000-10,000 µL	\$277.36	\$166.42

Special offers valid
January 19 – July 30, 2026