


Cell Culture Dishes



- Cell Culture Dishes are available in Tissue Culture Treated and Non-treated Cell Culture versions
- Petri (Non-Treated Cell Culture Dishes) have thicker and more durable plastic than standard Petri Dishes - both can be used for the same applications

- 1 60, 70 and 100mm dishes available with a pronounced grip ring for improved handling
- 2 Vented lids contain molded spacers to allow gas exchange (except Part No. 229655 which has a non-vented lid)

GRIP RING	PN 229660 & 229663	PN 229670	PN 229690, 229692 & 229693	PN 229620 & 229623
NO GRIP RING	PN 229635 & 229638	PN 229661	PN 229621	PN 229656
				PN 229650, 229651 & 229653

DISHES - TISSUE CULTURE TREATED

S, R/D F, NP, USP, PS


Part No.	Product Description	Grip Ring	Working Volume	Growth Area	Pack Description	Qty/Case
229635	35mm x 10mm Tissue Culture Treated Dish	N	5mL	8.5cm ²	10/Re-sealable Bag	500
229660	60mm x 15mm Tissue Culture Treated Dish	Y	7 - 8mL	21cm ²	10/Re-sealable Bag	500
NEW 229661	60mm x 15mm Tissue Culture Treated Dish	N	7 - 8mL	21cm ²	20/Re-sealable Bag	500
229670	70mm x 15mm Tissue Culture Treated Dish	Y	7 - 8mL	36.3cm ²	10/Re-sealable Bag	500
229690	100mm x 15mm Tissue Culture Treated Dish	Y	15 - 16mL	55cm ²	10/Re-sealable Bag	500
229620	100mm x 20mm Tissue Culture Treated Dish	Y	15 - 16mL	60.5cm ²	10/Re-sealable Bag	300
NEW 229621	100mm x 20mm Tissue Culture Treated Dish	N	15 - 16mL	57cm ²	10/Re-sealable Bag	300
229650	150mm x 20mm Tissue Culture Treated Dish	N	25 - 28mL	143cm ²	1/Bag	60
229651	150mm x 20mm Tissue Culture Treated Dish	N	25 - 28mL	143cm ²	5/Re-sealable Bag	100

DISHES - PETRI (NON-TREATED CELL CULTURE)


S, R/D F, NP, USP, PS

Part No.	Product Description	Grip Ring	Working Volume	Growth Area	Pack Description	Qty/Case
229638	35mm x 10mm Petri Dish	N	5mL	8.5cm ²	10/Bag	960
229663	60mm x 15mm Petri Dish	Y	7 - 8mL	21cm ²	10/Bag	500
229692	100mm x 15mm Petri Dish	Y	15 - 16mL	55cm ²	10/Bag	300
229693	100mm x 15mm Petri Dish	Y	15 - 16mL	55cm ²	10/Bag	500
229623	100mm x 20mm Petri Dish	Y	15 - 16mL	55cm ²	10/Bag	300
NEW 229656	150mm x 15mm Petri Dish	N	25 - 28mL	153cm ²	10/Bag	100
229653	150mm x 20mm Petri Dish	N	25 - 28mL	143cm ²	5/Bag	100
229655	150mm x 20mm Dish, Non-Vented Lid	N	25 - 28mL	143cm ²	5/Bag	100


Petri Dishes




Part No. 229665



Part No. 229695



Part No. 229694



Part No. 229697

- CELLTREAT Petri Dishes are manufactured in high capacity automated production lines resulting in high quality, value priced products
- Ventilation ribs reduce condensation and allow air circulation
- Consistent flatness allows for even media distribution

Petri Dish Selection Guide

STACKABLE	STACKABLE/SLIDEABLE	SLIPPABLE
Lids contain stacking rings for increased stability	Features partial edge ring on base which allows a single dish to be slid from a large stack and enhances stability when stacked	Designed for automated filling and streaking systems, features side arrows and ISO Mark Target for automation

DISHES - PETRI

S, R/D F, NP, USP, PS

Part No.	Product Description	Working Volume	Growth Area	Pack Description	Qty/Case	Dish Style
229665	60mm x 15mm Petri Dish	7 - 8 mL	21cm ²	20/Bag	500	Stackable
229694	100mm x 15mm Petri Dish, Slideable	15 - 16mL	55cm ²	20/Bag	500	Stackable/Slideable
229695	100mm x 15mm Petri Dish	15 - 16mL	55cm ²	20/Bag	500	Stackable
229697	100mm x 15mm Petri Dish, Slippable	15 - 16mL	55cm ²	20/Bag	500	Slippable

Compartment Petri Dishes

- Designed for simultaneous studies of organisms in a common environment
- Molded dividers create completely separate compartments while allowing airflow
- Ridge on bottom of dish allows easy handling and stacking
- Lids contain 3 molded spacers to allow gas exchange
- Working volume and growth area listed in table below are per compartment



DISHES - COMPARTMENT PETRI

S, R/D F, NP, USP, PS

Part No.	Product Description	Working Volume	Growth Area	Pack Description	Qty/Case	Stackable (Y/N)
229682	100mm x 15mm Petri Dish, 2 Compartments	7-8mL	28.3cm ²	20/Bag	500	Y
229683	100mm x 15mm Petri Dish, 3 Compartments	5mL	18.9cm ²	20/Bag	500	Y
229684	100mm x 15mm Petri Dish, 4 Compartments	4mL	14.1cm ²	20/Bag	500	Y