



- Tissue Culture Treated Flasks available in 6 sizes: 12.5, 25, 75, 182, 300 and **NEW** 225cm²
- Cell Culture Flasks are available in Tissue Culture Treated and Non-treated Suspension Culture versions
- Tissue Culture Treated Flasks have uniform hydrophilic surface ideal for consistent cell attachment
- Stable base enables filling with flask standing upright
- 182cm² flask offers more growth area than standard 175cm² flasks
- Convenient zip-closure bag is a practical alternative to folding and taping opened bags

- 1 Vent caps feature 0.22µm hydrophobic membrane
- 2 Large opening, canted neck design and moderate angles allow easy access to entire flask surface; ideal for cell scrapers (see page 21)
- 3 Molded, interlocking T design for improved stacking - included on 75, 182 & 300cm² and 250, 600 & 850mL size flasks only
- 4 Thick durable walls
- 5 Flask is printed with lot number and expiration date for easy identification and tracking

FLASKS - TISSUE CULTURE TREATED

S, R/D F, NP, USP, PS

Part No.	Product Description	Cap Style	Max Volume	Pack Description	Qty/Case
229320	12.5cm ² Tissue Culture Flask	Plug Seal	25mL	10/Re-sealable Bag	200
229330	25cm ² Tissue Culture Flask	Plug Seal	50mL	10/Re-sealable Bag	200
229340	75cm ² Tissue Culture Flask	Plug Seal	250mL	5/Re-sealable Bag	100
229350	182cm ² Tissue Culture Flask	Plug Seal	600mL	5/Re-sealable Bag	40
NEW 229370	225cm ² Tissue Culture Flask	Plug Seal	700mL	5/Re-sealable Bag	25
NEW 229360	300cm ² Tissue Culture Flask	Plug Seal	850mL	3/Re-sealable Bag	18
229321	12.5cm ² Tissue Culture Flask	Vent	25mL	10/Re-sealable Bag	200
229331	25cm ² Tissue Culture Flask	Vent	50mL	10/Re-sealable Bag	200
229341	75cm ² Tissue Culture Flask	Vent	250mL	5/Re-sealable Bag	100
229351	182cm ² Tissue Culture Flask	Vent	600mL	5/Re-sealable Bag	40
NEW 229371	225cm ² Tissue Culture Flask	Vent	700mL	5/Re-sealable Bag	25
229361	300cm ² Tissue Culture Flask	Vent	850mL	3/Re-sealable Bag	18



FLASKS - NON-TREATED SUSPENSION CULTURE

S, R/D F, NP, USP, PS

Part No.	Product Description	Cap Style	Max Volume	Pack Description	Qty/Case
229500	25mL Suspension Culture Flask	Vent	25mL	10/Re-sealable Bag	200
229510	50mL Suspension Culture Flask	Vent	50mL	10/Re-sealable Bag	200
229520	250mL Suspension Culture Flask	Vent	250mL	5/Re-sealable Bag	100
229530	600mL Suspension Culture Flask	Vent	600mL	5/Re-sealable Bag	40
NEW 229540	850mL Suspension Culture Flask	Vent	850mL	3/Re-sealable Bag	18

Flask Caps

FLASK CAPS

S, R/D F, NP, USP, HDPE

Part No.	Product Description	Cap Style	Fits Flask Size	Qty/Case
229392	Flask Cap	Plug Seal	12.5cm ² & 25mL	5
NEW 229389	Flask Cap	Plug Seal	25cm ² & 50mL	5
229396	Flask Cap	Plug Seal	75cm ² & 250mL	5
229398	Flask Cap	Plug Seal	182 & 300cm ² & 600 & 850mL	5
229393	Flask Cap	Vent	12.5cm ² & 25mL	5
NEW 229390	Flask Cap	Vent	25cm ² & 50mL	5
229397	Flask Cap	Vent	75cm ² & 250mL	5
229399	Flask Cap	Vent	182 & 300cm ² & 600 & 850mL	5