


## NEW Glass Bottom Cell Culture Dishes & Plates



- Glass bottom dishes and plates have the same great features as standard tissue culture treated products with the added benefit of borosilicate glass bottom microwells
- Convenient live cell viewing without needing to transfer
- Optically clear borosilicate glass for exceptional viewing
- Applications include: cell culture and imaging, live cell observation and specialty microscopy such as confocal high resolution and contrast
- Traditional footprints cause minimal disruption in cell culture techniques
- Constructed with noncytotoxic medical adhesive
- Round cover glass also available for use with glass bottom plates and dishes (see below)

DISHES - GLASS BOTTOM TISSUE CULTURE TREATED		S, R/D F, NP, USP, PS, GLS	
Part No.	Product Description	Pack Description	Qty/Case
229632	30mm x 10mm Tissue Culture Treated Dish, 15mm Glass Bottom	10/Resealable Bag	50
PLATES - GLASS BOTTOM TISSUE CULTURE TREATED		S, R/D F, NP, USP, PS, GLS	
Part No.	Product Description	Pack Description	Qty/Case
229107	6 Well Tissue Culture Plate with Lid, 20mm Glass Bottom	Individual	5
229125	24 Well Tissue Culture Plate with Lid, 10mm Glass Bottom	Individual	5


## Round Cover Glass

- High-quality, optically clear borosilicate glass cover slips are ideal for specialty microscopy and imaging
- Ideal for use with Dishes and Multiple Well Plates and anywhere a square cover glass is unsuitable
- Consistent flatness and scratch-free
- Resistant to organic solvents
- Packaged in a hinged plastic storage case
- 15mm Round Cover Glass can be used with a 24 Well Glass Bottom Plate (Part No. 229125)
- 20mm Round Cover Glass can be used with a 12 Well Plate or 30mm Glass Bottom Dish (Part No. 229632)
- 25mm Round Cover Glass can be used with a 6 Well Glass Bottom Plate (Part No. 229107)



ROUND COVER GLASS		S, R/D F, NP, USP, GLS	
Part No.	Product Description	Qty/Case	
229172	15mm Round Cover Glass, Fits 24 Well Plate	100/Pack	
229173	20mm Round Cover Glass, Fits 12 Well Plate or 30mm Dish	100/Pack	
229174	25mm Round Cover Glass, Fits 6 Well Plate	100/Pack	

## NEW Chambered Cell Culture Slides



- Grow cells directly on the slide without any need for future transfer; ideal for growing and viewing cells under a microscope on standard 25mm x 75mm slide
- Uniform hydrophilic surface ideal for consistent cell attachment
- Slides are labeled with numeric markings
- Packing trays can be used as incubation racks in CO<sub>2</sub> incubators
- Adhesive free

- 1 Integrated gasket prevents cross-contamination between chambers
- 2 Inert hydrophobic green border printed on slide to define cell culture area
- 3 Chamber holder features easy flip locks for assembly/disassembly; no separation tool required

CHAMBERED CELL CULTURE SLIDES		S, R/D F, NP, USP, PS, PP, GLS	
Part No.	Product Description	Pack Description	Qty/Case
229161	1 Chamber Cell Culture Slide, Glass	6/Tray with 2 Trays/Box	12
229162	2 Chamber Cell Culture Slide, Glass	6/Tray with 2 Trays/Box	12
229164	4 Chamber Cell Culture Slide, Glass	6/Tray with 2 Trays/Box	12
229168	8 Chamber Cell Culture Slide, Glass	6/Tray with 2 Trays/Box	12