

Erlenmeyer & Fernbach Flasks



- Designed for mammalian suspension culture applications as well as bacterial & yeast culture
- Conveniently pre-sterilized
- 1 70mm opening for easy filling and dispensing (Erlenmeyer & Fernbach)
- 2 Over-sized vent area in cap for improved gas exchange (Erlenmeyer & Fernbach)
- 3 Fernbach flask shape increases surface area to volume ratio
- 4 Baffled bottom improves mixing and aeration
 - Plain bottom also available

Part No. 229896 - Transfer Assembly

ERLENMEYER & FERNBACH FLASKS S, NP, USP, A, PC

Part No.	Product Description	Cap Type	Bottom Type	Max. Volume	Pack Description	Qty/Case
229850	2L Erlenmeyer Flask	Vent	Plain	2200mL	1/Bag	6
229855	2L Erlenmeyer Flask	Vent	Baffled	2200mL	1/Bag	6
229860	2L Fernbach Flask	Vent	Plain	2400mL	1/Bag	6
229865	2L Fernbach Flask	Vent	Baffled	2400mL	1/Bag	6
229870	3L Fernbach Flask	Vent	Plain	4000mL	1/Bag	4
229875	3L Fernbach Flask	Vent	Baffled	4000mL	1/Bag	4


ERLENMEYER AND FERNBACH FLASK CAPS S, R/D F, NP, USP, HDPE

Part No.	Product Description	Cap Type	Bottom Type	Max. Volume	Pack Description	Qty/Case
229894	CAP ONLY, 70mm, Erlenmeyer/Fernbach	Solid	NA	NA	1/Bag	24
229895	CAP ONLY, 70mm, Erlenmeyer/Fernbach	Vent	NA	NA	1/Bag	24

TRANSFER ASSEMBLY S

NEW 229896	Transfer Assembly for 3L Fernbach Flask	NA	NA	NA	1/Bag	1
-------------------	---	----	----	----	-------	---

Solution Bottles



- Compatible with bottle top filters or for use as solution storage containers (Bottle Top Filters sold separately, see page 35)

- 1 Solution bottle neck size is 45mm
- 2 Curved bottle design for easy handling

SOLUTION BOTTLES		S, R/D F, NP, USP, PS		
Part No.	Product Description	Working volume	Pack Description	Qty/Case
229781	150mL Solution Bottle	150mL	1/Bag	24
229782	250mL Solution Bottle	250mL	1/Bag	24
229784	500mL Solution Bottle	500mL	1/Bag	24
229785	1000mL Solution Bottle	1000mL	1/Bag	24

NEW 500mL Square Media Bottle

- Ideal for media prep and storage of laboratory reagents and biological buffers
- Square design saves valuable shelf space and is ideal for stacking, such as in a refrigerator
- Unique individually bagged packaging allows opening one at a time
- PETG material reduces gas permeability

- 1 Molded graduations on two sides
- 2 Flat surface polyethylene (HDPE) caps are ideal for marking

SQUARE MEDIA BOTTLE S, R/D F, NP, USP, PETG

Part No.	Product Description	Working volume	Pack Description	Qty/Case
229788	500mL Media Bottle, Square	500mL	1/Bag	24