

Centrifuge Bottles

- Centrifuge Bottles have a rigid body for low to moderate speed centrifugation of biological and chemical samples
- High quality polypropylene has excellent chemical resistance
- Centrifuge Bottles are reusable and autoclavable
- Bottle must be filled 75% of total capacity to obtain maximum results
- Packaged non-sterile, 2 per re-sealable bag




CENTRIFUGE TUBES		N, USP, A, PP	
Part No.	Product Description	Qty/Case	Max Force (RCF) g
229467	250mL Centrifuge Bottle	2	27500
229468	500mL Centrifuge Bottle	2	13500

Non-Treated Roller Bottles

- Non-Treated Roller Bottles can be used for suspension cell culture and large capacity liquid storage
- 1 Large knurled cap for easy grip
- 2 Easy-to-read black printed graduations



Vent caps available separately for use with open systems (Part No. 229391)



NON-TREATED ROLLER BOTTLES		S, R/D F, NP, USP, PS		
Part No.	Product Description	Cap Style	Pack Description	Qty/Case
229582	1000mL Roller Bottle	Non-Vented	1/Bag	24
229584	2000mL Roller Bottle	Non-Vented	1/Bag	12
ROLLER BOTTLES - CAP ONLY		S, R/D F, NP, USP, HDPE		
229391	Roller Bottle, Cap Only	Vented	1/Bag	24

CELLTREAT Scientific Products

The CELLTREAT Scientific Products team understands that quality, innovation, and value are important to our customers. The entire organization is dedicated to providing quality cell culture and liquid handling products, exceptional customer service, and measurable value to the scientific research community.

We believe that innovation is a critical component to providing exceptional customer service. The CELLTREAT product offering includes many products with unique and differentiated features that can only be found under the CELLTREAT brand name. We strive to bring new products to market quickly without sacrificing quality at a price that represents a real value for our customers.

www.celltreat.com



Visit us on the web and save time!

We make it easy for you to access all of the information you need in one central location:

- Expanded product information including the latest new products
- Request a Quality Certificate (COC)
- Request product samples
- Current CELLTREAT news and promotions
- Contact us

Quality Information



CELLTREAT Scientific Products can provide a comprehensive Quality Certificate for all products. The CELLTREAT Quality Certificate includes the Part Number, Lot Number, Sterility Assurance Level, Endotoxin/Pyrogen Free Assurance, Non-Cytotoxicity Assurance and RNase/DNase Free Assurance. The CELLTREAT Quality Certificate was designed to include those items found on Certificates of Compliance and Certificates of Assurance. The Quality Certificate is specific to each part number and lot of manufactured product.

To request a copy of a Quality Certificate, please complete the online request form at: www.celltreat.com/quality-certificate-request

Samples

Need to request a sample?

Simply send an email to samples@celltreat.com or complete the online sample request form at www.celltreat.com/sample-request, include your name and complete address along with the part number of the product you are requesting. We will send you an email with tracking information as soon as it has shipped. It's that easy!

