

## Centrifuge Tubes - Bag



- Convenient zip-closure bag is a practical alternative to folding and taping opened bags
- Available in 15 & 50mL Conical Bottom, 50mL Self-Standing, 15 & 50mL Flip Top and 15 & 50mL Black tube options
- Easily refill reusable racks with tubes packed in bags

### CENTRIFUGE TUBES - BAG

S, R/D F, NP, USP, A, PP

Part No.	Product Description	Pack Description	Qty/Case	Sterile (Y/N)	Max Force (RCF) g
229411	15mL Centrifuge Tube	25/Re-sealable Bag	500	Y	12000
229421	50mL Centrifuge Tube	25/Re-sealable Bag	500	Y	16000
229430	50mL Centrifuge Tube - Self-Standing	25/Re-sealable Bag	500	Y	6000

### FLIP TOP CENTRIFUGE TUBES - BAG

S, R/D F, NP, USP, A, PP

Part No.	Product Description	Pack Description	Qty/Case	Sterile (Y/N)	Max Force (RCF) g
229417	15mL Centrifuge Tube, Flip Top Cap	25/Re-sealable Bag	500	Y	9400
229427	50mL Centrifuge Tube, Flip Top Cap	25/Re-sealable Bag	500	Y	9400
229425	50mL Centrifuge Tube, Flip Top Cap Only	25/Re-sealable Bag	1000	Y	NA

### BLACK CENTRIFUGE TUBES - BAG

S, R/D F, NP, USP, A, PP


Part No.	Product Description	Pack Description	Qty/Case	Sterile (Y/N)	Max Force (RCF) g
229431	15mL Centrifuge Tube, Black	25/Re-sealable Bag	100	Y	12000
229432	15mL Centrifuge Tube, Black	25/Re-sealable Bag	500	Y	12000
229433	50mL Centrifuge Tube, Black	25/Re-sealable Bag	100	Y	16000
229434	50mL Centrifuge Tube, Black	25/Re-sealable Bag	500	Y	16000

### CAP OR TUBE ONLY CENTRIFUGE TUBES - BAG

N, R/D F, NP, USP, PP

Part No.	Product Description	Pack Description	Qty/Case	Sterile (Y/N)	Max Force (RCF) g
229450	TUBE ONLY, 15mL Centrifuge Tube	500/Bag	500	N	12000
229451	CAP ONLY, 15mL Centrifuge Tube	500/Bag	500	N	NA
229452	TUBE & CAP, 15mL Centrifuge Tube & Cap	Tubes 500/Bag & Caps 500/Bag	500	N	12000
229453	CAP ONLY, 50mL Centrifuge Tube	500/Bag	500	N	NA
229454	TUBE ONLY, 50mL Centrifuge Tube	500/Bag	500	N	16000
229456	TUBE & CAP, 50mL Centrifuge Tube & Cap	Tubes 500/Bag & Caps 500/Bag	500	N	16000

## Centrifuge Tubes - Foam Rack




- Sturdy alternative to disposable paperboard racks
- Also available with printed marking spot and pre-numbered caps

### CENTRIFUGE TUBES - FOAM RACK

S, R/D F, NP, A, PP

Part No.	Product Description	Pack Description	Qty/Case	Sterile (Y/N)	Max Force (RCF) g
229412	15mL Centrifuge Tube	50/Foam Rack	300	Y	12000
229422	50mL Centrifuge Tube	25/Foam Rack	300	Y	16000
229423	50mL Centrifuge Tube, Marking Spot Cap	25/Foam Rack	300	Y	16000
229424	50mL Centrifuge Tube, Numbered Cap	25/Foam Rack	300	Y	16000

## Centrifuge Tubes - Plastic Rack



▲ Reusable/Recyclable Racks also available for purchase separately (PN 229419 & 229429)


- Plastic Racks are manufactured from polypropylene and are autoclavable, reusable and recyclable
- Chemically resistant
- Design of the tray prevents racks from sticking/freezing to surfaces
- Widely spaced design catches drips and allows tubes to be retrieved more easily

### CENTRIFUGE TUBES - PLASTIC RACK

R/D F, NP, USP, A, PP

Part No.	Product Description	Pack Description	Qty/Case	Sterile (Y/N)	Max Force (RCF) g
229418	15mL Centrifuge Tube	25/Plastic Rack with 2 Racks/Bag	300	Y	12000
229428	50mL Centrifuge Tube	25/Plastic Rack	300	Y	16000
229419	PLASTIC RACK ONLY, 15mL Centrifuge Tube	5/Bag	5	N	NA
229429	PLASTIC RACK ONLY, 50mL Centrifuge Tube	5/Bag	5	N	NA

## Centrifuge Tubes - Paperboard Rack



- Recyclable alternative to disposable foam rack
- Available in 15 & 50mL Conical Bottom and 15 & 50mL Flip Top tube options
- Folds flat for disposal or recycling

### CENTRIFUGE TUBES - PAPERBOARD RACK

R/D F, NP, USP, A, PP

Part No.	Product Description	Pack Description	Qty/Case	Sterile (Y/N)	Max Force (RCF) g
229410	15mL Centrifuge Tube	25/Paperboard Rack with 2 Racks/Bag	500	Y	12000
229420	50mL Centrifuge Tube	25/Paperboard Rack	500	Y	16000
229455	50mL Centrifuge Tube	25/Paperboard Rack	500	N	16000

### FLIP TOP CENTRIFUGE TUBES - PAPERBOARD RACK

S, R/D F, NP, USP, A, PP

Part No.	Product Description	Pack Description	Qty/Case	Sterile (Y/N)	Max Force (RCF) g
229416	15mL Centrifuge Tube, Flip Top Cap	25/Paperboard Rack with 2 Racks/Bag	500	Y	9400
229426	50mL Centrifuge Tube, Flip Top Cap	25/Paperboard Rack	500	Y	9400