



#### **Rotor FA-45-24-11**

- Aerosol-Tight
- Max. 24 × 1.5/2.0 mL tubes
- Max. 21,130 × g (15,000 rpm)



#### **Rotor FA-45-24-11-Special**

- Aerosol-Tight
- Max. 24 × 1.5/2.0 mL tubes
- Max. 21,130 × g (15,000 rpm)
- PTFE coating for increased chemical resistance



#### **Rotor F-45-18-11-Kit**

- Max. 18 × spin columns or 1.5/2.0 mL tubes
- Max. 18,111 × g (15,000 rpm)
- Extended rim that protects tube caps from shearing off



#### **Rotor F-45-32-5-PCR**

- Max. 4 × PCR strips
- Max. 18,615 × g (15,000 rpm)