

hybex™

Media Storage Bottles



**Borosilicate
Glass**

hybex™ reusable media bottles are manufactured from high quality, borosilicate 3.3 glass. With their superior chemical resistance, these bottles are ideal for storage of reagents, culture media, biological fluids and a variety of other aqueous and nonaqueous solutions.

A low coefficient of thermal expansion (3.3), provides temperature resistance that is far superior to that of standard lab glass (or soda lime glass). Hybex glass can withstand temperatures from -70 to 500°C, permitting repeated use in freezers, microwaves, autoclaves and dry heat sterilization chambers.

All hybex bottles have permanent, white enamel graduations for volume measurement and a white identification box for marking or coding.

All hybex media bottles include a drip-free sealing ring as well as a standard blue polypropylene (GL45) sealing cap. Bottles are also available with green, purple, white and yellow caps, or in an assortment pack containing two of each color.

- ✦ Borosilicate 3.3 hybex glass
- ✦ Autoclavable caps (GL45) & drip rings
- ✦ Permanent enamel graduation markings
- ✦ Large area for bottle marking/identification
- ✦ Autoclavable, reuseable & breakage resistant
- ✦ Five cap colors

Technical Data

Material:	hybex™ borosilicate glass (3.3)
Cap Thread Size:	GL45*
Temp. Range:	-70°C to 500°C
Sterilization:	Autoclave or dry heat
Cap Material:	Polypropylene (140°C)
Cap Colors:	Blue, green, purple yellow, white or assorted

*50ml bottles are only supplied with blue GL32 sealing caps

Ordering Information

B3000-50	Media storage bottle, 50ml with standard (GL32) blue cap, 10/pk.
B3000-100*	Media storage bottle, 100ml with standard (GL45) blue cap, 10/pk.
B3000-250*	Media storage bottle, 250ml with standard (GL45) blue cap, 10/pk.
B3000-500*	Media storage bottle, 500ml with standard (GL45) blue cap, 10/pk.
B3000-1000*	Media storage bottle, 1L with standard (GL45) blue cap, 10/pk.
B3000-2000*	Media storage bottle, 2L with standard (GL45) blue cap, 5/pk.
B3000-SP*	Media storage bottle, Starter Pack with standard (GL45) blue cap (Starter Pack of ten bottles includes 2x100ml, 3x250ml, 3x500ml & 2x1000ml)
B3000-CAP*	High temperature red caps, up to 180°C (GL45), 10/pk
B3000-HTC	Replacement cap, standard blue (GL45), 10/pk.
B3000-RIN	Replacement sealing ring, (GL45), 10/pk.
B3000-CAP2	Replacement cap, standard blue (GL32), 10/pk.
B3000-RIN2	Replacement sealing ring, (GL32), 10/pk.

*Please add (-G) for green caps, (-P) for purple caps, (-Y) for yellow caps, (-W) for white caps or (-AS) for assorted caps (2 of ea. color)

Benchmark
Scientific



Benchmark Scientific Inc.
PO Box 709
Edison, NJ 08818
Phone: 908-769-5555
Fax: 908-222-1864

Web: www.BenchmarkScientific.com
Email: Info@BenchmarkScientific.com