


**NEW! Glass Bottom Cell Culture Dishes and Plates**



- Glass bottom dishes and plates have the same great features as standard tissue culture treated products with the added benefit of borosilicate glass bottom microwells
- Convenient live cell viewing without needing to transfer
- Optically clear borosilicate glass for exceptional viewing
- Applications include: cell culture and imaging, live cell observation and specialty microscopy such as confocal high resolution and contrast
- Constructed with noncytotoxic medical adhesive
- Round cover glass also available for use with glass bottom plates and dishes

Part No.	Product Description	Qty/Case
<b>Glass Bottom Tissue Culture Dish</b>		
229632	30mm x 10mm TCT Dish, 15mm Glass Bottom	50
<b>Glass Bottom Tissue Culture Well Plates</b>		
229107	6 Well TCT Plate with Lid, 20mm Glass Bottom	5
229125	24 Well TCT Plate with Lid, 10mm Glass Bottom	5



**NEW! Round Cover Glass**

- High-quality, optically clear borosilicate glass cover slips are ideal for specialty microscopy and imaging
- Thickness: #1.5 (0.16mm-0.19mm)
- Sterile



Part No.	Product Description	Compatible With	Qty/Case
229172	15mm Round Cover Glass, Fits 24 Well Plate	Standard 24 Multiple Well Plates and 24 Well Glass Bottom Plate (229125)	100
229173	20mm Round Cover Glass, Fits 12 Well Plate or 30mm Dish	Standard 12 Multiple Well Plates and 30mm Glass Bottom Dish (229632)	100
229174	25mm Round Cover Glass, Fits 6 Well Plate	Standard 6 Multiple Well and Glass Bottom Plate (229107)	100

**NEW! Chambered Cell Culture Slides**



- Grow cells directly on the slide without any need for future transfer; ideal for growing and viewing cells under a microscope on standard 25mm x 75mm slide
- Uniform hydrophilic surface ideal for consistent cell attachment
- Slides are labeled with numeric markings
- Packing trays can be used as incubation racks in CO<sub>2</sub> incubators
- Adhesive free

- 1 Integrated gasket prevents cross-contamination between chambers
- 2 Inert hydrophobic green border printed on slide to define cell culture area
- 3 Chamber holder features easy flip locks for assembly/disassembly; no separation tool required

CELLTREAT Part No.	Product Description	Pack Description	Qty/Case
229161	1 Chamber Cell Culture Slide, Glass	6/Tray with 2 Trays/Box	12
229162	2 Chamber Cell Culture Slide, Glass	6/Tray with 2 Trays/Box	12
229164	4 Chamber Cell Culture Slide, Glass	6/Trays with 2 Trays/Box	12
229168	8 Chamber Cell Culture Slide, Glass	6/Trays with 2 Trays/Box	12