


# NEW Best Value Serological Pipets



SEROLOGICAL PIPETS - INDIVIDUALLY WRAPPED (PAPER/PLASTIC) IN BAGS				S, R/D F, NP, USP, PS	
Part No.	Product Description	Pack Description	Qty/Case	Negative Graduations	
229001B	1mL Serological Pipet	100 Individually Wrapped/Bag	500	-0.3mL	
229002B	2mL Serological Pipet	100 Individually Wrapped/Bag	500	-0.5mL	
229005B	5mL Serological Pipet	50 Individually Wrapped/Bag	200	-3mL	
229010B	10mL Serological Pipet	50 Individually Wrapped/Bag	200	-3mL	
229212B	10mL Serological Pipet, Wide Tip	50 Individually Wrapped/Bag	200	-3mL	
229025B	25mL Serological Pipet	50 Individually Wrapped/Bag	200	-10mL	
229030B	50mL Serological Pipet	25 Individually Wrapped/Bag	100	-10mL	

SEROLOGICAL PIPETS - BULK PACKED IN BAGS				S, R/D F, NP, USP, PS	
Part No.	Product Description	Pack Description	Qty/Case	Negative Graduations	
229231	1mL Serological Pipet	25 Bulk Packed/Re-sealable Bag	1000	-0.3mL	
229232	2mL Serological Pipet	25 Bulk Packed/Re-sealable Bag	800	-0.5mL	
229235	5mL Serological Pipet	25 Bulk Packed/Re-sealable Bag	500	-3mL	
229240	10mL Serological Pipet	25 Bulk Packed/Re-sealable Bag	500	-3mL	
229245	25mL Serological Pipet	25 Bulk Packed/Re-sealable Bag	200	-10mL	
229250	50mL Serological Pipet	25 Bulk Packed/Re-sealable Bag	200	-10mL	

# NEW Best Value Serological Pipets



10mL Standard vs. 10mL Wide Tip Pipet

All serological pipets are manufactured to the ASTM Standard Specification for Serological Pipets E-934

- 1 Bi-directional (positive and negative) graduations are printed in bold easy-to-read black ink
- 2 10mL Wide Tip Serological Pipet for viscous solutions - Wide Tip opening is 4x larger than standard 10mL pipet
- 3 25 & 50mL pipets feature an anti-drip tip for accurate delivery
- 4 Premium material barrier plug
- 5 Color coded stripes and packaging enable quick volume identification
- 6 Printed indicator arrow for quickly locating negative graduations when using a pipet controller

1, 2, 5 & 10mL Serological Pipets



## Serological Pipets - Plastic/Plastic Wrap

- Peel-to-open plastic/plastic wrappers are the ideal option for eliminating particulate; proven benefit over traditional paper/plastic wrapper
- Wrappers reduce static cling versus other brands of plastic wrap pipets
- All plastic wrapper can be wiped with alcohol without smearing ink or penetrating paper
- Same great features as standard pipets



### PLASTIC/PLASTIC WRAP SEROLOGICAL PIPETS - INDIVIDUALLY WRAPPED IN BAGS

S, R/D F, NP, USP, PS

Part No.	Product Description	Pack Description	Qty/Case
229201A	1mL Serological Pipet	100 Individually Wrapped/Bag	800
229202A	2mL Serological Pipet	100 Individually Wrapped/Bag	600
229205A	5mL Serological Pipet	50 Individually Wrapped/Bag	200
229210A	10mL Serological Pipet	50 Individually Wrapped/Bag	200
229225A	25mL Serological Pipet	50 Individually Wrapped/Bag	200
229230A	50mL Serological Pipet	30 Individually Wrapped/Bag	90

## Serological Pipets - Short



Short vs. standard length pipet comparison

- Compact design of short serological pipet provides ergonomic benefit and reduces storage space
- Shorter length helps prevent touching pipet tip to surfaces in the hood
- Ideal for students and those new to pipetting
- Same great features as standard length pipets

### SEROLOGICAL PIPETS - SHORT

S, R/D F, NP, USP, PS

Part No.	Product Description	Pack Description	Qty/Case
229206B	5mL Serological Pipet, Short	50 Individually Wrapped/Bag	200
229211B	10mL Serological Pipet, Short	50 Individually Wrapped/Bag	200
229226B	25mL Serological Pipet, Short	25 Individually Wrapped/Bag	100

## Serological Pipets - Bacteriological/Milk

- Bacteriological/Milk pipets are manufactured to meet APHA requirements shown in Standard Methods for Examination of Dairy Products and ASTM E1043-00(2010) Standard Specification for Pipets, Milk and Cream Examination
- Specialized graduations allow quick volume identification
- Same great features as standard pipets

### SEROLOGICAL PIPETS - BACTERIOLOGICAL/MILK

S, R/D F, NP, USP, PS

Part No.	Product Description	Pack Description	Qty/Case
229267	1.1mL Bacteriological/Milk Pipet	50 Individually Wrapped/Bag	250
229257	1.1mL Bacteriological/Milk Pipet	50 Bulk Packed/Bag	500
229268	2.2mL Bacteriological/Milk Pipet	50 Individually Wrapped/Bag	250
229258	2.2mL Bacteriological/Milk Pipet	50 Bulk Packed/Bag	400



# Classic Serological Pipets



- 1 Color coded stripes and packaging enable quick volume identification
- 2 Bi-directional graduations are printed in bold easy-to-read black ink
- 3 Premium material filter plug

### Assorted packaging options:

- 4 Individually wrapped in bags
- 5 Individually wrapped in convenient shelf carton
- 6 Individually wrapped in plastic/plastic wrapper (see page 14 for details)

## SEROLOGICAL PIPETS - INDIVIDUALLY WRAPPED IN BAGS S, R/D F, NP, USP, PS



Part No.	Product Description	Pack Description	Qty/Case
229201B	1mL Serological Pipet	100 Individually Wrapped/Bag	800
229202B	2mL Serological Pipet	100 Individually Wrapped/Bag	600
229205B	5mL Serological Pipet	50 Individually Wrapped/Bag	200
229210B	10mL Serological Pipet	50 Individually Wrapped/Bag	200
229225B	25mL Serological Pipet	50 Individually Wrapped/Bag	200
229230B	50mL Serological Pipet	30 Individually Wrapped/Bag	90
229220B	100mL Serological Pipet	10 Individually Wrapped/Bag	50

## SEROLOGICAL PIPETS - INDIVIDUALLY WRAPPED IN DISPENSER CARTON S, R/D F, NP, USP, PS



Part No.	Product Description	Pack Description	Qty/Case
229201	1mL Serological Pipet	200 Individually Wrapped/Dispenser Carton	800
229202	2mL Serological Pipet	150 Individually Wrapped/Dispenser Carton	600
229205	5mL Serological Pipet	100 Individually Wrapped/Dispenser Carton	200
229210	10mL Serological Pipet	100 Individually Wrapped/Dispenser Carton	200
229225	25mL Serological Pipet	50 Individually Wrapped/Dispenser Carton	200
229230	50mL Serological Pipet	45 Individually Wrapped/Dispenser Carton	90

# Serological Pipets - Open End



- Ideal for viscous and high particulate solutions
- Wide, open end is ideal for slurries such as Sephadex
- Same great features as standard length pipets

## SEROLOGICAL PIPETS - OPEN END S, R/D F, NP, USP, PS

Part No.	Product Description	Pack Description	Qty/Case
229223B	5mL Serological Pipet, Open End	50 Individually Wrapped/Bag	200
229243	5mL Serological Pipet, Open End	25 Bulk Packed/Bag	500
229224B	10mL Serological Pipet, Open End	50 Individually Wrapped/Bag	200
229244	10mL Serological Pipet, Open End	25 Bulk Packed/Bag	500