


## NEW CELLTREAT Freeze Cryogenic Vial Storage Boxes

- Durable polycarbonate material designed to withstand a wide temperature range and is a reusable, autoclavable alternative to disposable racks
- Ventilation holes in the lids and bases prevent vacuum sealing in freezer
- Storage Boxes are labeled with numeric markings for easy identification of individual positions
- Safe to use at ultra-low temperatures
- 100 place rack features ridges along the sides to improve handling



### CELLTREAT FREEZE CRYOGENIC VIAL STORAGE BOXES

N, R/D F, NP, USP, A, PC

Part No.	Product Description	Compatible With	Qty/Case
229942	Storage Box, CF Cryogenic Vial, 25 Place	Fits 1mL-1.8mL vials	5
229943	Storage Box, CF Cryogenic Vial, 81 Place	Fits 1mL-1.8mL vials	5
229945	Storage Box, Cryogenic Vial, 100 Place	Fits 1mL-1.8mL vials	5
229944	Deep Storage Box, CF Cryogenic Vial, 81 Place	Fits 4.5mL vials	5

## NEW CELLTREAT Freeze Cryogenic Vials Cap Inserts



- Easily color-code samples with colored cap inserts
- CELLTREAT Freeze Cap Inserts are compatible with all CELLTREAT Freeze Cryogenic Vials

### CELLTREAT FREEZE CRYOGENIC VIALS - CAP INSERTS

N, R/D F, NP, USP, PP

Part No.	Product Description	Color	Pack Description	Qty/Case
229930	Cap Insert for CF Cryogenic Vials	Red	500 Bulk Packed/Re-sealable Bag	500
229931	Cap Insert for CF Cryogenic Vials	Orange	500 Bulk Packed/Re-sealable Bag	500
229932	Cap Insert for CF Cryogenic Vials	Yellow	500 Bulk Packed/Re-sealable Bag	500
229933	Cap Insert for CF Cryogenic Vials	Green	500 Bulk Packed/Re-sealable Bag	500
229934	Cap Insert for CF Cryogenic Vials	Blue	500 Bulk Packed/Re-sealable Bag	500
229935	Cap Insert for CF Cryogenic Vials	Purple	500 Bulk Packed/Re-sealable Bag	500
229936	Cap Insert for CF Cryogenic Vials	Pink	500 Bulk Packed/Re-sealable Bag	500
229937	Cap Insert for CF Cryogenic Vials	Brown	500 Bulk Packed/Re-sealable Bag	500
229938	Cap Insert for CF Cryogenic Vials	Gray	500 Bulk Packed/Re-sealable Bag	500
229939	Cap Insert for CF Cryogenic Vials	White	500 Bulk Packed/Re-sealable Bag	500
229940	Cap Insert for CF Cryogenic Vials	Assorted	500 Bulk Packed/Re-sealable Bag	500

## Roller Bottles



- Available in two sizes: 1000mL volume with 490cm<sup>2</sup> of surface area and 2000mL volume with 850cm<sup>2</sup> of surface area, treated or non-treated versions
- Expanded Surface Roller Bottles (ESRB) also available
- Vented or non-vented cap for open or closed systems
- Tissue Culture Treated Roller Bottles have uniform hydrophilic surface ideal for consistent cell attachment

- 1 Large knurled cap for easy grip
- 2 Ribbed surface of ESRB increases surface area from 850cm<sup>2</sup> to 1900cm<sup>2</sup>
- 3 Easy-to-read black printed graduations



Part No. 229387 - 1900cm<sup>2</sup> ESRB, Vent Cap

### ROLLER BOTTLES - TISSUE CULTURE TREATED

S, R/D F, NP, USP, PS

Part No.	Product Description	Cap Style	Pack Description	Qty/Case
229382	490cm <sup>2</sup> Roller Bottle	Non-Vented	1/Bag	24
229383	490cm <sup>2</sup> Roller Bottle	Vented	1/Bag	24
229384	850cm <sup>2</sup> Roller Bottle	Non-Vented	1/Bag	12
229385	850cm <sup>2</sup> Roller Bottle	Vented	1/Bag	12

### ROLLER BOTTLES - TISSUE CULTURE TREATED, EXPANDED SURFACE

S, R/D F, NP, USP, PS

229386	1900cm <sup>2</sup> Expanded Surface Roller Bottle	Non-Vented	1/Bag	12
229387	1900cm <sup>2</sup> Expanded Surface Roller Bottle	Vented	1/Bag	12

### ROLLER BOTTLES - NON-TREATED

S, R/D F, NP, USP, PS

229582	1000mL Roller Bottle	Non-Vented	1/Bag	24
229584	2000mL Roller Bottle	Non-Vented	1/Bag	12

### ROLLER BOTTLES - CAP ONLY

S, R/D F, NP, USP, HDPE

229391	Roller Bottle, Cap Only	Vented	1/Bag	24
--------	-------------------------	--------	-------	----