

## Markers

- Fast drying, smudge resistant, extra-fine point marker
- Ideal for permanently marking all plastic lab ware, especially centrifuge tubes



### MARKERS

Part No.	Product Description	Qty/Case
229405	Green Tube Marker, Fast Drying Ink	5
229406	Blue Tube Marker, Fast Drying Ink	5
229407	Black Tube Marker, Fast Drying Ink	5
229408	Red Tube Marker, Fast Drying Ink	5

## CELLTREAT Accessories



### BAG CUTTER

- Safely and easily cuts plastic bags and pouches
- Features unique contoured design, sharp non-removable stainless steel blade and lanyard hang hole for convenience
- Includes integrated tape splitter for easy box opening

### TIMER


- Multiple functions include timer, clock and alarm
- Includes easel stand, mounting magnet and retractable hook for hanging
- Uses 1 standard AAA Battery (not included)

### ACCESSORIES

Part No.	Product Description	Qty/Case
230100	Bag Cutter	5
230123	Multi-Function Timer	1

## Loops, Needles and Spreaders

- Flexible for easy access to small volume tubes and dishes
- Color coded for easy identification
- Gentle on delicate agar or other surfaces
- L-Shaped Cell Spreader is ideal for evenly spreading bacterial cultures
- Loops and L-Shaped Cell Spreader available individually wrapped in peel-to-open paper/plastic wrappers



### L-SHAPED CELL SPREADERS

S, R/D F, NP, USP, PP

Part No.	Product Description	Color	Length	Pack Description	Qty/Case
229617	L-Shaped Spreader	Green	145mm	10/Pack	500
229616	L-Shaped Spreader	Green	145mm	100 Individually Wrapped/Bag	100

### INOCULATING LOOPS AND NEEDLES

S, R/D F, NP, USP, PP

229601	1uL Inoculating Loop	White	218mm	20/Pack	2000
229602	1uL Inoculating Loop	White	218mm	100 Individually Wrapped/Bag	600
229610	10uL Inoculating Loop	Blue	221mm	20/Pack	2000
229612	10uL Inoculating Loop	Blue	221mm	100 Individually Wrapped/Bag	600
229615	Inoculating Needle	Yellow	221mm	20/Pack	2000

## Pipet Controller

- Lightweight (6.3oz/170g) and designed to fit comfortably in either hand
- Long-life rechargeable lithium battery provides up to 8 hours of use on a two hour charge


- 1 Adjustable dial for speed setting
- 2 Cordless and rechargeable
- 3 Low pressure fingertip control buttons
- 4 Easy fit adapter accepts pipet sizes from 1 to 100mL



### PIPET CONTROLLER

Part No.	Product Description	Pack Description	Qty/Case
230200	Pipet Controller, 120V	Individual	1
230202	AC Adapter, Replacement, 120V	Individual	1
230204	Battery, Replacement	Individual	1
230205	Rubber Grommet, Replacement	Individual	1
230206	Filter, Replacement, 0.45µm	2/Bag	2

## Pipet Storage Rack



- Durable ABS plastic rack can stand on flat surfaces as well as be attached to the metal sides of lab cabinets and appliances

- 1 Strong magnets are included on both sides for increased versatility
- 2 Angled shelves provide secure yet easy access to pipets
- 3 Pipet sizes, colors and types can be viewed from the side windows and front
- 4 Side bracket can hold pipet fillers, bulbs and pumps

### PIPET STORAGE RACK

ABS

Part No.	Product Description	Qty/Case
230295	Pipet Storage Rack	1

## Reagent Reservoirs

- Reagent Reservoirs are designed to assist in repetitive liquid transfer
- White color enhances solution visibility
- Angled design promotes funneling of liquid to maximize retrieval



### REAGENT RESERVOIRS

S, R/D F, NP, USP, PS

Part No.	Product Description	Pack Description	Color	Qty/Case
229290	50mL Reagent Reservoir	1 Individually Wrapped/Bag	White	80
229291	50mL Reagent Reservoir	5 Bulk Packed/Bag	White	200