

# CAPP *Maestro*

ELECTRONIC PIPETTES

**CAPP *Maestro*** is the electronic pipette that combines all benefits that the modern researcher needs. The Favourites button is the first function that you will see on **CAPP *Maestro***'s fully coloured OLED display. The coloured OLED display gives better optical control while pipetting as it has better contrast and a greater screen refresh rate compared to the other LCD screens. Also, an OLED display consumes less power compared to the other displays.

Lightweight construction and elegant ergonomic design with a wide hook allows you to rest the pipette at the palm of your hand with no need to actually "grab it".

**CAPP *Maestro***'s tip ejector also works as a finger rest for your thumb. Simply rest your thumb on the tip ejector, and gently activate the pipetting button at the same time.

Optimize your serial work with advanced functions such as multi dispensing, dilution and mixing, all available on **CAPP *Maestro*** while easily navigating through the user-friendly graphic intuitive interface that is self explanatory and can be operated very easily!

For your routine lab work, **CAPP *Maestro*** allows you to easily store your favorite protocols and to switch from one to the other in a matter of 2 clicks.

If you are tired of exchanging the batteries too often, then try **CAPP *Maestro*** that uses lithium batteries which are ideal for long-lasting performance. You can also use this electronic pipette while it is charging.

Calibration and password setting functions are available on this pipette.



## Ordering Information

Cat. No.	Description	Inaccuracy %	Imprecision %
M10-1	CAPPMaestro single channel, 0,2-10 $\mu$ L	2.50 / 1.00	1.80 / 0.50
M100-1	CAPPMaestro single channel, 5-100 $\mu$ L	2.50 / 0.70	1.00 / 0.20
M200-1	CAPPMaestro single channel, 10-200 $\mu$ L	2.50 / 0.60	0.70 / 0.20
M1000-1	CAPPMaestro single channel, 50-1000 $\mu$ L	2.00 / 0.80	0.60 / 0.20
M5000-1	CAPPMaestro single channel, 500-5000 $\mu$ L	2.40 / 0.60	0.60 / 0.15