



- ① Adjusting mechanism (PP)
- ② Piston seat (PP)
- ③ Outer sleeve (PP)
- ④ Volume selector (PP, PEEK)
- ⑤ Piston (PFA, ETFE)
- ⑥ Cylinder sleeve (PP)
- ⑦ Glass cylinder (Borosilicate glass)
- ⑧ Valve head (PFA), housing(PP)
- ⑨ Valve ball (Borosilicate glass)
- ⑩ Intake valve (ETFE/PTFE)
- ⑪ Telescopic intake tube (FEP/ETFE)
- ⑫ Plug (PP)
- ⑬ Valve switch (PP)
- ⑭ Discharge tube support (PP)
- ⑮ Discharge/recirculation valve (PFA/PTFE)
- ⑯ Discharge tube securing nut (PP)
- ⑰ Discharge tube (FEP)
- ⑱ Discharge tube closure cap (PVDF)
- ⑲ Tool (Size2.5-10ml, size25-100ml)
- ⑳ Stop ring(PP)
- ㉑ Securing clip