

## COMPATIBILITY CHART

| Description   | Cat. No. | CAPP   |                       |               | Gilson          |                        |                      | Eppendorf   |  |                       | Brand         |                 |                        | Thermo               |             |  | Sartorius             |               |                 | HTL                    |                      |             |           |             |              |            |               |   |
|---|----------|--|-----------------------|---------------|-----------------|------------------------|----------------------|-------------|--|-----------------------|---------------|-----------------|------------------------|----------------------|-------------|--|-----------------------|---------------|-----------------|------------------------|----------------------|-------------|-----------|-------------|--------------|------------|---------------|---|
|   |          | Bravo / ecopipette / CAPPaero / CAPP Maestro / CAPP Ironic | Pipetman / Pipetman G | Research Plus | Transferpette S | Fimmpipette F1 / Novus | Proline Plus / Picus | Labmate Pro | Bravo / ecopipette / CAPPaero / CAPP Maestro / CAPP Ironic | Pipetman / Pipetman G | Research Plus | Transferpette S | Fimmpipette F1 / Novus | Proline Plus / Picus | Labmate Pro | Bravo / ecopipette / CAPPaero / CAPP Maestro / CAPP Ironic | Pipetman / Pipetman G | Research Plus | Transferpette S | Fimmpipette F1 / Novus | Proline Plus / Picus | Labmate Pro |           |             |              |            |               |   |
| <b>CAPP ExpellPlus Low Retention Pipette Tips</b>               |          | 0.2-2 µl   | 0.5-10 µl             | 20-200 µl     | 100-1000 µl     | 1000-10000 µl          | 0.2-2 µl             | 0.5-10 µl   | 20-200 µl  | 100-1000 µl           | 1000-10000 µl | 0.5-10 µl pink  | 2-20 µl turquoise      | 2-20 µl yellow       | 100-1000 µl | 0.5-10 µl  | 20-200 µl             | 100-1000 µl   | 1000-10000 µl   | 0.2-10 µl              | 0.5-10 µl            | 10-100 µl   | 20-200 µl | 100-1000 µl | 1000-5000 µl | 50-1000 µl | 1000-10000 µl |   |
| ExpellPlus 10 µL, clear, sterile                                | 5030010C | ■  | ■                     | ■             | ■               | ■                      | ■                    | ■           | ■  | ■                     | ■             | ■               | ■                      | ■                    | ■           | ■  | ■                     | ■             | ■               | ■                      | ■                    | ■           | ■         | ■           | ■            | ■          | ■             |   |
| ExpellPlus 10 µL, clear, sterile                                | 5030020C | ■  | ■                     | ■             | ■               | ■                      | ■                    | ■           | ■  | ■                     | ■             | ■               | ■                      | ■                    | ■           | ■  | ■                     | ■             | ■               | ■                      | ■                    | ■           | ■         | ■           | ■            | ■          | ■             |   |
| ExpellPlus 10 µL XL, clear, sterile                             | 5030040C | ■  | ■                     | ■             | ■               | ■                      | ■                    | ■           | ■  | ■                     | ■             | ■               | ■                      | ■                    | ■           | ■  | ■                     | ■             | ■               | ■                      | ■                    | ■           | ■         | ■           | ■            | ■          | ■             |   |
| ExpellPlus 10 µL XL, clear, sterile                             | 5030050C | ■  | ■                     | ■             | ■               | ■                      | ■                    | ■           | ■  | ■                     | ■             | ■               | ■                      | ■                    | ■           | ■  | ■                     | ■             | ■               | ■                      | ■                    | ■           | ■         | ■           | ■            | ■          | ■             |   |
| ExpellPlus 10/20 µL XL, extra narrow, clear, sterile            | 5030051C | ■  | ■                     | ■             | ■               | ■                      | ■                    | ■           | ■  | ■                     | ■             | ■               | ■                      | ■                    | ■           | ■  | ■                     | ■             | ■               | ■                      | ■                    | ■           | ■         | ■           | ■            | ■          | ■             |   |
| ExpellPlus 200 µL, clear, sterile                               | 5030070C | ■  | ■                     | ■             | ■               | ■                      | ■                    | ■           | ■  | ■                     | ■             | ■               | ■                      | ■                    | ■           | ■  | ■                     | ■             | ■               | ■                      | ■                    | ■           | ■         | ■           | ■            | ■          | ■             |   |
| ExpellPlus 200 µL, clear, sterile                               | 5030080C | ■  | ■                     | ■             | ■               | ■                      | ■                    | ■           | ■  | ■                     | ■             | ■               | ■                      | ■                    | ■           | ■  | ■                     | ■             | ■               | ■                      | ■                    | ■           | ■         | ■           | ■            | ■          | ■             |   |
| ExpellPlus 300 µL, clear, sterile                               | 5130110C | ■  | ■                     | ■             | ■               | ■                      | ■                    | ■           | ■  | ■                     | ■             | ■               | ■                      | ■                    | ■           | ■  | ■                     | ■             | ■               | ■                      | ■                    | ■           | ■         | ■           | ■            | ■          | ■             |   |
| <b>CAPP ExpellPlus Low Retention Filter Tips</b>                |          | 0.2-2 µl   | 0.5-10 µl             | 20-200 µl     | 100-1000 µl     | 1000-10000 µl          | 0.2-2 µl             | 0.5-10 µl   | 20-200 µl  | 100-1000 µl           | 1000-10000 µl | 0.5-10 µl pink  | 2-20 µl turquoise      | 2-20 µl yellow       | 100-1000 µl | 0.5-10 µl  | 20-200 µl             | 100-1000 µl   | 1000-10000 µl   | 0.2-10 µl              | 0.5-10 µl            | 10-100 µl   | 20-200 µl | 100-1000 µl | 1000-5000 µl | 50-1000 µl | 1000-10000 µl |   |
| ExpellPlus 10 µL, w/ filter, clear, sterile                     | 5030030C | ■  | ■                     | ■             | ■               | ■                      | ■                    | ■           | ■  | ■                     | ■             | ■               | ■                      | ■                    | ■           | ■  | ■                     | ■             | ■               | ■                      | ■                    | ■           | ■         | ■           | ■            | ■          | ■             |   |
| ExpellPlus 10 µL XL, w/ filter, clear, sterile                  | 5030060C | ■  | ■                     | ■             | ■               | ■                      | ■                    | ■           | ■  | ■                     | ■             | ■               | ■                      | ■                    | ■           | ■  | ■                     | ■             | ■               | ■                      | ■                    | ■           | ■         | ■           | ■            | ■          | ■             |   |
| ExpellPlus 10/20 µL XL, extra narrow, w/ filter, clear, sterile | 5030061C | ■  | ■                     | ■             | ■               | ■                      | ■                    | ■           | ■  | ■                     | ■             | ■               | ■                      | ■                    | ■           | ■  | ■                     | ■             | ■               | ■                      | ■                    | ■           | ■         | ■           | ■            | ■          | ■             |   |
| ExpellPlus 20 µL, w/ filter, clear, sterile                     | 5030062C | ■  | ■                     | ■             | ■               | ■                      | ■                    | ■           | ■  | ■                     | ■             | ■               | ■                      | ■                    | ■           | ■  | ■                     | ■             | ■               | ■                      | ■                    | ■           | ■         | ■           | ■            | ■          | ■             |   |
| ExpellPlus 100 µL, w/ filter, clear, sterile                    | 5030066C | ■  | ■                     | ■             | ■               | ■                      | ■                    | ■           | ■  | ■                     | ■             | ■               | ■                      | ■                    | ■           | ■  | ■                     | ■             | ■               | ■                      | ■                    | ■           | ■         | ■           | ■            | ■          | ■             |   |
| ExpellPlus 200 µL, w/ filter, clear, sterile                    | 5030090C | ■  | ■                     | ■             | ■               | ■                      | ■                    | ■           | ■  | ■                     | ■             | ■               | ■                      | ■                    | ■           | ■  | ■                     | ■             | ■               | ■                      | ■                    | ■           | ■         | ■           | ■            | ■          | ■             |   |
| ExpellPlus 300 µL, w/ filter, clear, sterile                    | 5030115C | ■  | ■                     | ■             | ■               | ■                      | ■                    | ■           | ■  | ■                     | ■             | ■               | ■                      | ■                    | ■           | ■  | ■                     | ■             | ■               | ■                      | ■                    | ■           | ■         | ■           | ■            | ■          | ■             |   |
| ExpellPlus 1000 µL, w/ filter, clear, sterile                   | 5030150C | ■  | ■                     | ■             | ■               | ■                      | ■                    | ■           | ■  | ■                     | ■             | ■               | ■                      | ■                    | ■           | ■  | ■                     | ■             | ■               | ■                      | ■                    | ■           | ■         | ■           | ■            | ■          | ■             |   |
| <b>CAPP Expell Pipette Tips</b>                                 |          | 0.2-2 µl   | 0.5-10 µl             | 20-200 µl     | 100-1000 µl     | 1000-10000 µl          | 0.2-2 µl             | 0.5-10 µl   | 20-200 µl  | 100-1000 µl           | 1000-10000 µl | 0.5-10 µl pink  | 2-20 µl turquoise      | 2-20 µl yellow       | 100-1000 µl | 0.5-10 µl  | 20-200 µl             | 100-1000 µl   | 1000-10000 µl   | 0.2-10 µl              | 0.5-10 µl            | 10-100 µl   | 20-200 µl | 100-1000 µl | 1000-5000 µl | 50-1000 µl | 1000-10000 µl |   |
| Expell 10 µL, clear   | 5030015C | ■  | ■                     | ■             | ■               | ■                      | ■                    | ■           | ■  | ■                     | ■             | ■               | ■                      | ■                    | ■           | ■  | ■                     | ■             | ■               | ■                      | ■                    | ■           | ■         | ■           | ■            | ■          | ■             | ■ |
| Expell 10 µL, clear   | 5130010C | ■  | ■                     | ■             | ■               | ■                      | ■                    | ■           | ■  | ■                     | ■             | ■               | ■                      | ■                    | ■           | ■  | ■                     | ■             | ■               | ■                      | ■                    | ■           | ■         | ■           | ■            | ■          | ■             | ■ |
| Expell 10 µL XL, clear  | 5030045C | ■  | ■                     | ■             | ■               | ■                      | ■                    | ■           | ■  | ■                     | ■             | ■               | ■                      | ■                    | ■           | ■  | ■                     | ■             | ■               | ■                      | ■                    | ■           | ■         | ■           | ■            | ■          | ■             | ■ |
| Expell 10 µL XL, clear  | 5130040C | ■  | ■                     | ■             | ■               | ■                      | ■                    | ■           | ■  | ■                     | ■             | ■               | ■                      | ■                    | ■           | ■  | ■                     | ■             | ■               | ■                      | ■                    | ■           | ■         | ■           | ■            | ■          | ■             | ■ |
| Expell 10/20 µL XL, extra narrow, clear                         | 5030046C | ■  | ■                     | ■             | ■               | ■                      | ■                    | ■           | ■  | ■                     | ■             | ■               | ■                      | ■                    | ■           | ■  | ■                     | ■             | ■               | ■                      | ■                    | ■           | ■         | ■           | ■            | ■          | ■             | ■ |
| Expell 10/20 µL XL, extra narrow, clear                         | 5130041C | ■  | ■                     | ■             | ■               | ■                      | ■                    | ■           | ■  | ■                     | ■             | ■               | ■                      | ■                    | ■           | ■  | ■                     | ■             | ■               | ■                      | ■                    | ■           | ■         | ■           | ■            | ■          | ■             | ■ |
| Expell 200 µL, clear  | 5030075C | ■  | ■                     | ■             | ■               | ■                      | ■                    | ■           | ■  | ■                     | ■             | ■               | ■                      | ■                    | ■           | ■  | ■                     | ■             | ■               | ■                      | ■                    | ■           | ■         | ■           | ■            | ■          | ■             | ■ |
| Expell 200 µL, clear  | 5130070C | ■  | ■                     | ■             | ■               | ■                      | ■                    | ■           | ■  | ■                     | ■             | ■               | ■                      | ■                    | ■           | ■  | ■                     | ■             | ■               | ■                      | ■                    | ■           | ■         | ■           | ■            | ■          | ■             | ■ |
| Expell 200 µL, yellow   | 4130075C | ■  | ■                     | ■             | ■               | ■                      | ■                    | ■           | ■  | ■                     | ■             | ■               | ■                      | ■                    | ■           | ■  | ■                     | ■             | ■               | ■                      | ■                    | ■           | ■         | ■           | ■            | ■          | ■             | ■ |
| Expell 300 µL, clear  | 5130105C | ■  | ■                     | ■             | ■               | ■                      | ■                    | ■           | ■  | ■                     | ■             | ■               | ■                      | ■                    | ■           | ■  | ■                     | ■             | ■               | ■                      | ■                    | ■           | ■         | ■           | ■            | ■          | ■             | ■ |
| Expell 300 µL, clear  | 5130100C | ■  | ■                     | ■             | ■               | ■                      | ■                    | ■           | ■  | ■                     | ■             | ■               | ■                      | ■                    | ■           | ■  | ■                     | ■             | ■               | ■                      | ■                    | ■           | ■         | ■           | ■            | ■          | ■             | ■ |
| Expell 1000 µL (1250 µL), clear                                 | 5130135C | ■  | ■                     | ■             | ■               | ■                      | ■                    | ■           | ■  | ■                     | ■             | ■               | ■                      | ■                    | ■           | ■  | ■                     | ■             | ■               | ■                      | ■                    | ■           | ■         | ■           | ■            | ■          | ■             | ■ |
| Expell 1000 µL (1250 µL), clear                                 | 5130130C | ■  | ■                     | ■             | ■               | ■                      | ■                    | ■           | ■  | ■                     | ■             | ■               | ■                      | ■                    | ■           | ■  | ■                     | ■             | ■               | ■                      | ■                    | ■           | ■         | ■           | ■            | ■          | ■             | ■ |
| Expell 1000 µL (1250 µL), blue                                  | 4130135C | ■  | ■                     | ■             | ■               | ■                      | ■                    | ■           | ■  | ■                     | ■             | ■               | ■                      | ■                    | ■           | ■  | ■                     | ■             | ■               | ■                      | ■                    | ■           | ■         | ■           | ■            | ■          | ■             | ■ |
| Expell 1000 µL (1250 µL), clear, sterile                        | 5130140C | ■  | ■                     | ■             | ■               | ■                      | ■                    | ■           | ■  | ■                     | ■             | ■               | ■                      | ■                    | ■           | ■  | ■                     | ■             | ■               | ■                      | ■                    | ■           | ■         | ■           | ■            | ■          | ■             | ■ |
| Expell 1250 µL XL, clear  | 5130121C | ■  | ■                     | ■             | ■               | ■                      | ■                    | ■           | ■  | ■                     | ■             | ■               | ■                      | ■                    | ■           | ■  | ■                     | ■             | ■               | ■                      | ■                    | ■           | ■         | ■           | ■            | ■          | ■             | ■ |
| Expell 1250 µL XL, clear  | 5130120C | ■  | ■                     | ■             | ■               | ■                      | ■                    | ■           | ■  | ■                     | ■             | ■               | ■                      | ■                    | ■           | ■  | ■                     | ■             | ■               | ■                      | ■                    | ■           | ■         | ■           | ■            | ■          | ■             | ■ |
| Expell 1250 µL XL, clear, sterile                               | 5130122C | ■  | ■                     | ■             | ■               | ■                      | ■                    | ■           | ■  | ■                     | ■             | ■               | ■                      | ■                    | ■           | ■  | ■                     | ■             | ■               | ■                      | ■                    | ■           | ■         | ■           | ■            | ■          | ■             | ■ |
| Expell 5000 µL, clear, sterile                                  | 5130165C | ■  | ■                     | ■             | ■               | ■                      | ■                    | ■           | ■  | ■                     | ■             | ■               | ■                      | ■                    | ■           | ■  | ■                     | ■             | ■               | ■                      | ■                    | ■           | ■         | ■           | ■            | ■          | ■             | ■ |
| Expell 5000 µL, clear   | 5130160C | ■  | ■                     | ■             | ■               | ■                      | ■                    | ■           | ■  | ■                     | ■             | ■               | ■                      | ■                    | ■           | ■  | ■                     | ■             | ■               | ■                      | ■                    | ■           | ■         | ■           | ■            | ■          | ■             | ■ |
| Expell 10000 µL, clear, sterile                                 | 5130195C | ■  | ■                     | ■             | ■               | ■                      | ■                    | ■           | ■  | ■                     | ■             | ■               | ■                      | ■                    | ■           | ■  | ■                     | ■             | ■               | ■                      | ■                    | ■           | ■         | ■           | ■            | ■          | ■             | ■ |
| Expell 10000 µL, clear  | 5130190C | ■  | ■                     | ■             | ■               | ■                      | ■                    | ■           | ■  | ■                     | ■             | ■               | ■                      | ■                    | ■           | ■  | ■                     | ■             | ■               | ■                      | ■                    | ■           | ■         | ■           | ■            | ■          | ■             | ■ |
| Expell 10000 µL, clear  | 5130191C | ■  | ■                     | ■             | ■               | ■                      | ■                    | ■           | ■  | ■                     | ■             | ■               | ■                      | ■                    | ■           | ■  | ■                     | ■             | ■               | ■                      | ■                    | ■           | ■         | ■           | ■            | ■          | ■             | ■ |

This compatibility chart is for your reference. We are not responsible for technical modifications or changes made to the mentioned pipettes other than those of CAPP. Specifications like inaccuracy and precision values have not been checked. For a complete compatibility list please contact [info@cappahn.com](mailto:info@cappahn.com).

## COMPATIBILITY CHART

| Description   | Cat. No. | CAPP  |                       |               |                   | Gilson                 |                      |             |   | Eppendorf             |               |                   |                        | Brand                |             |   |                       | Thermo        |                   |                        |                      | Sartorius   |   |                       |               | HTL               |                        |                      |             |   |
|---|----------|---|-----------------------|---------------|-------------------|------------------------|----------------------|-------------|---|-----------------------|---------------|-------------------|------------------------|----------------------|-------------|---|-----------------------|---------------|-------------------|------------------------|----------------------|-------------|---|-----------------------|---------------|-------------------|------------------------|----------------------|-------------|---|
|   |          | Bravo / ecopipette / CAPP Aero / CAPP Maestro / CAPP Tronic | Pipetman / Pipetman G | Research Plus | Transferpipette S | Finnpipette E1 / Novus | Proline Plus / Picus | Labmate Pro | Bravo / ecopipette / CAPP Aero / CAPP Maestro / CAPP Tronic | Pipetman / Pipetman G | Research Plus | Transferpipette S | Finnpipette E1 / Novus | Proline Plus / Picus | Labmate Pro | Bravo / ecopipette / CAPP Aero / CAPP Maestro / CAPP Tronic | Pipetman / Pipetman G | Research Plus | Transferpipette S | Finnpipette E1 / Novus | Proline Plus / Picus | Labmate Pro | Bravo / ecopipette / CAPP Aero / CAPP Maestro / CAPP Tronic | Pipetman / Pipetman G | Research Plus | Transferpipette S | Finnpipette E1 / Novus | Proline Plus / Picus | Labmate Pro |   |
| <b>CAPP Expell Filter Tips</b>                                      |          |   |                       |               |                   |                        |                      |             |   |                       |               |                   |                        |                      |             |   |                       |               |                   |                        |                      |             |   |                       |               |                   |                        |                      |             |   |
| Expell 10 µL, w/ filter, clear, sterile                             | 5130030C | ■   | ■                     | ■             | ■                 | ■                      | ■                    | ■           | ■   | ■                     | ■             | ■                 | ■                      | ■                    | ■           | ■   | ■                     | ■             | ■                 | ■                      | ■                    | ■           | ■   | ■                     | ■             | ■                 | ■                      | ■                    | ■           |   |
| Expell 10 µL, w/ filter, clear                                      | 5131030C |   | ■                     | ■             | ■                 | ■                      | ■                    | ■           | ■   | ■                     | ■             | ■                 | ■                      | ■                    | ■           | ■   | ■                     | ■             | ■                 | ■                      | ■                    | ■           | ■   | ■                     | ■             | ■                 | ■                      | ■                    | ■           | ■ |
| Expell 10 µL XL, w/ filter, clear, sterile                          | 5130060C | ■   | ■                     | ■             | ■                 | ■                      | ■                    | ■           | ■   | ■                     | ■             | ■                 | ■                      | ■                    | ■           | ■   | ■                     | ■             | ■                 | ■                      | ■                    | ■           | ■   | ■                     | ■             | ■                 | ■                      | ■                    | ■           | ■ |
| Expell 10 µL XL, w/ filter, clear                                   | 5131060C |   | ■                     | ■             | ■                 | ■                      | ■                    | ■           | ■   | ■                     | ■             | ■                 | ■                      | ■                    | ■           | ■   | ■                     | ■             | ■                 | ■                      | ■                    | ■           | ■   | ■                     | ■             | ■                 | ■                      | ■                    | ■           | ■ |
| Expell 10/20 µL XL, extra narrow, w/ filter, clear, sterile         | 5130061C | ■   | ■                     | ■             | ■                 | ■                      | ■                    | ■           | ■   | ■                     | ■             | ■                 | ■                      | ■                    | ■           | ■   | ■                     | ■             | ■                 | ■                      | ■                    | ■           | ■   | ■                     | ■             | ■                 | ■                      | ■                    | ■           | ■ |
| Expell 10/20 µL XL, extra narrow, w/ filter, clear                  | 5131061C |   | ■                     | ■             | ■                 | ■                      | ■                    | ■           | ■   | ■                     | ■             | ■                 | ■                      | ■                    | ■           | ■   | ■                     | ■             | ■                 | ■                      | ■                    | ■           | ■   | ■                     | ■             | ■                 | ■                      | ■                    | ■           | ■ |
| Expell 20 µL, w/ filter, clear, sterile                             | 5130062C | ■   | ■                     | ■             | ■                 | ■                      | ■                    | ■           | ■   | ■                     | ■             | ■                 | ■                      | ■                    | ■           | ■   | ■                     | ■             | ■                 | ■                      | ■                    | ■           | ■   | ■                     | ■             | ■                 | ■                      | ■                    | ■           | ■ |
| Expell 20 µL, w/ filter, clear                                      | 5131062C |   | ■                     | ■             | ■                 | ■                      | ■                    | ■           | ■   | ■                     | ■             | ■                 | ■                      | ■                    | ■           | ■   | ■                     | ■             | ■                 | ■                      | ■                    | ■           | ■   | ■                     | ■             | ■                 | ■                      | ■                    | ■           | ■ |
| Expell 100 µL, w/ filter, clear, sterile                            | 5130066C | ■   | ■                     | ■             | ■                 | ■                      | ■                    | ■           | ■   | ■                     | ■             | ■                 | ■                      | ■                    | ■           | ■   | ■                     | ■             | ■                 | ■                      | ■                    | ■           | ■   | ■                     | ■             | ■                 | ■                      | ■                    | ■           | ■ |
| Expell 100 µL, w/ filter, clear                                     | 5131066C |   | ■                     | ■             | ■                 | ■                      | ■                    | ■           | ■   | ■                     | ■             | ■                 | ■                      | ■                    | ■           | ■   | ■                     | ■             | ■                 | ■                      | ■                    | ■           | ■   | ■                     | ■             | ■                 | ■                      | ■                    | ■           | ■ |
| Expell 200 µL, w/ filter, clear, sterile                            | 5130090C | ■   | ■                     | ■             | ■                 | ■                      | ■                    | ■           | ■   | ■                     | ■             | ■                 | ■                      | ■                    | ■           | ■   | ■                     | ■             | ■                 | ■                      | ■                    | ■           | ■   | ■                     | ■             | ■                 | ■                      | ■                    | ■           | ■ |
| Expell 200 µL, w/ filter, clear                                     | 5131090C |   | ■                     | ■             | ■                 | ■                      | ■                    | ■           | ■   | ■                     | ■             | ■                 | ■                      | ■                    | ■           | ■   | ■                     | ■             | ■                 | ■                      | ■                    | ■           | ■   | ■                     | ■             | ■                 | ■                      | ■                    | ■           | ■ |
| Expell 1000 µL, w/ filter, clear, sterile                           | 5130150C | ■   | ■                     | ■             | ■                 | ■                      | ■                    | ■           | ■   | ■                     | ■             | ■                 | ■                      | ■                    | ■           | ■   | ■                     | ■             | ■                 | ■                      | ■                    | ■           | ■   | ■                     | ■             | ■                 | ■                      | ■                    | ■           | ■ |
| Expell 1000 µL, w/ filter, clear                                    | 5131150C |   | ■                     | ■             | ■                 | ■                      | ■                    | ■           | ■   | ■                     | ■             | ■                 | ■                      | ■                    | ■           | ■   | ■                     | ■             | ■                 | ■                      | ■                    | ■           | ■   | ■                     | ■             | ■                 | ■                      | ■                    | ■           | ■ |
| Expell 1250 µL XL, w/ filter, clear, sterile                        | 5130123C | ■   | ■                     | ■             | ■                 | ■                      | ■                    | ■           | ■   | ■                     | ■             | ■                 | ■                      | ■                    | ■           | ■   | ■                     | ■             | ■                 | ■                      | ■                    | ■           | ■   | ■                     | ■             | ■                 | ■                      | ■                    | ■           | ■ |
| <b>CAPP Expell and ExpellPlus individually wrapped Pipette Tips</b> |          |   |                       |               |                   |                        |                      |             |   |                       |               |                   |                        |                      |             |   |                       |               |                   |                        |                      |             |   |                       |               |                   |                        |                      |             |   |
| ExpellPlus 10 µL, clear, sterile, indiv. wrapped                    | 5030012C | ■   | ■                     | ■             | ■                 | ■                      | ■                    | ■           | ■   | ■                     | ■             | ■                 | ■                      | ■                    | ■           | ■   | ■                     | ■             | ■                 | ■                      | ■                    | ■           | ■   | ■                     | ■             | ■                 | ■                      | ■                    | ■           | ■ |
| ExpellPlus 10 µL XL, clear, sterile, indiv. wrapped                 | 5030042C |   | ■                     | ■             | ■                 | ■                      | ■                    | ■           | ■   | ■                     | ■             | ■                 | ■                      | ■                    | ■           | ■   | ■                     | ■             | ■                 | ■                      | ■                    | ■           | ■   | ■                     | ■             | ■                 | ■                      | ■                    | ■           | ■ |
| ExpellPlus 200 µL, clear, sterile, indiv. wrapped                   | 5030072C | ■   | ■                     | ■             | ■                 | ■                      | ■                    | ■           | ■   | ■                     | ■             | ■                 | ■                      | ■                    | ■           | ■   | ■                     | ■             | ■                 | ■                      | ■                    | ■           | ■   | ■                     | ■             | ■                 | ■                      | ■                    | ■           | ■ |
| Expell 1000 µL (1250 µL), clear, sterile, indiv. wrapped            | 5130132C | ■   | ■                     | ■             | ■                 | ■                      | ■                    | ■           | ■   | ■                     | ■             | ■                 | ■                      | ■                    | ■           | ■   | ■                     | ■             | ■                 | ■                      | ■                    | ■           | ■   | ■                     | ■             | ■                 | ■                      | ■                    | ■           | ■ |
| <b>CAPP Expell and ExpellPlus Filter Tips in PaperBox Packs</b>     |          |   |                       |               |                   |                        |                      |             |   |                       |               |                   |                        |                      |             |   |                       |               |                   |                        |                      |             |   |                       |               |                   |                        |                      |             |   |
| ExpellPlus 10 µL, w/ filter, clear, sterile                         | 5030029C | ■   | ■                     | ■             | ■                 | ■                      | ■                    | ■           | ■   | ■                     | ■             | ■                 | ■                      | ■                    | ■           | ■   | ■                     | ■             | ■                 | ■                      | ■                    | ■           | ■   | ■                     | ■             | ■                 | ■                      | ■                    | ■           | ■ |
| ExpellPlus 10/20 µL XL, extra narrow, w/ filter, clear, sterile     | 5030059C | ■   | ■                     | ■             | ■                 | ■                      | ■                    | ■           | ■   | ■                     | ■             | ■                 | ■                      | ■                    | ■           | ■   | ■                     | ■             | ■                 | ■                      | ■                    | ■           | ■   | ■                     | ■             | ■                 | ■                      | ■                    | ■           | ■ |
| ExpellPlus 20 µL, w/ filter, clear, sterile                         | 5030063C | ■   | ■                     | ■             | ■                 | ■                      | ■                    | ■           | ■   | ■                     | ■             | ■                 | ■                      | ■                    | ■           | ■   | ■                     | ■             | ■                 | ■                      | ■                    | ■           | ■   | ■                     | ■             | ■                 | ■                      | ■                    | ■           | ■ |
| ExpellPlus 100 µL, w/ filter, clear, sterile                        | 5030064C | ■   | ■                     | ■             | ■                 | ■                      | ■                    | ■           | ■   | ■                     | ■             | ■                 | ■                      | ■                    | ■           | ■   | ■                     | ■             | ■                 | ■                      | ■                    | ■           | ■   | ■                     | ■             | ■                 | ■                      | ■                    | ■           | ■ |
| ExpellPlus 200 µL, w/ filter, clear, sterile                        | 5030089C | ■   | ■                     | ■             | ■                 | ■                      | ■                    | ■           | ■   | ■                     | ■             | ■                 | ■                      | ■                    | ■           | ■   | ■                     | ■             | ■                 | ■                      | ■                    | ■           | ■   | ■                     | ■             | ■                 | ■                      | ■                    | ■           | ■ |
| ExpellPlus 300 µL, w/ filter, clear, sterile                        | 5030114C | ■   | ■                     | ■             | ■                 | ■                      | ■                    | ■           | ■   | ■                     | ■             | ■                 | ■                      | ■                    | ■           | ■   | ■                     | ■             | ■                 | ■                      | ■                    | ■           | ■   | ■                     | ■             | ■                 | ■                      | ■                    | ■           | ■ |
| Expell 1000 µL, w/ filter, clear, sterile                           | 5130149C | ■   | ■                     | ■             | ■                 | ■                      | ■                    | ■           | ■   | ■                     | ■             | ■                 | ■                      | ■                    | ■           | ■   | ■                     | ■             | ■                 | ■                      | ■                    | ■           | ■   | ■                     | ■             | ■                 | ■                      | ■                    | ■           | ■ |
| ExpellPlus 1250 µL XL, w/ filter, clear, sterile                    | 5130119C | ■   | ■                     | ■             | ■                 | ■                      | ■                    | ■           | ■   | ■                     | ■             | ■                 | ■                      | ■                    | ■           | ■   | ■                     | ■             | ■                 | ■                      | ■                    | ■           | ■   | ■                     | ■             | ■                 | ■                      | ■                    | ■           | ■ |
| <b>CAPP Expell Pipette Tips Reload System</b>                       |          |   |                       |               |                   |                        |                      |             |   |                       |               |                   |                        |                      |             |   |                       |               |                   |                        |                      |             |   |                       |               |                   |                        |                      |             |   |
| Expell 10 µL, clear, reload   | 5030035C | ■   | ■                     | ■             | ■                 | ■                      | ■                    | ■           | ■   | ■                     | ■             | ■                 | ■                      | ■                    | ■           | ■   | ■                     | ■             | ■                 | ■                      | ■                    | ■           | ■   | ■                     | ■             | ■                 | ■                      | ■                    | ■           | ■ |
| Expell 10 µL XL, extra narrow, clear, reload                        | 5030055C | ■   | ■                     | ■             | ■                 | ■                      | ■                    | ■           | ■   | ■                     | ■             | ■                 | ■                      | ■                    | ■           | ■   | ■                     | ■             | ■                 | ■                      | ■                    | ■           | ■   | ■                     | ■             | ■                 | ■                      | ■                    | ■           | ■ |
| Expell 200 µL, clear, reload  | 5030095C | ■   | ■                     | ■             | ■                 | ■                      | ■                    | ■           | ■   | ■                     | ■             | ■                 | ■                      | ■                    | ■           | ■   | ■                     | ■             | ■                 | ■                      | ■                    | ■           | ■   | ■                     | ■             | ■                 | ■                      | ■                    | ■           | ■ |
| Expell 300 µL, clear, reload  | 5130111C | ■   | ■                     | ■             | ■                 | ■                      | ■                    | ■           | ■   | ■                     | ■             | ■                 | ■                      | ■                    | ■           | ■   | ■                     | ■             | ■                 | ■                      | ■                    | ■           | ■   | ■                     | ■             | ■                 | ■                      | ■                    | ■           | ■ |
| Expell 1000 µL (1250 µL), clear, reload                             | 5130155C | ■   | ■                     | ■             | ■                 | ■                      | ■                    | ■           | ■   | ■                     | ■             | ■                 | ■                      | ■                    | ■           | ■   | ■                     | ■             | ■                 | ■                      | ■                    | ■           | ■   | ■                     | ■             | ■                 | ■                      | ■                    | ■           | ■ |

This compatibility chart is for your reference. We are not responsible for technical modifications or changes made to the mentioned pipettes other than those of CAPP. Specifications like inaccuracy and imprecision values have not been checked. For a complete compatibility list please contact [info@cappahn.com](mailto:info@cappahn.com).

■ compatible

■ not compatible

