

Oxford LTS Tips LTS/XLS Pipette Tips

Features and Benefits

Advantages

- Most cost-effective LTS tip on the market
- Designed to fit LTS/XLS Rainin pipettes
- Ultra fine point tip ends to guarantee accuracy, precision, and consistency

Substance Free

- RNase, DNase, cytotoxic, PCR inhibitors, endotoxic, pyrogenic free
- Substance-free hydrophobic resins to prevent contamination and sample retention

High Quality Manufacturing

- Manufactured with diamond finished molds using high quality Polypropylene material
- Built to fit requiring minimal force to attach and detach pipettes.
- QC tested and regulated to adhere to Oxford high quality standards
- Low retention, Sterilized, Filtered



Product Specifications for Oxford Lab Products LTS Tips

Catalog Numbers	LTR-20-SLF, LTR-200-SLF, LTR-1000-SLF
Specification	Non-Pyrogenic, DNase & RNase, Nucleic Acid, Trace Metal, PCR inhibitor, and ATP free
Qty / Case	96 tips/rack, 10 racks/pack, 5 packs/case
Material	High Quality Polypropylene and Hydrophobic Resins
Sterile	Yes
Low Retention	Yes