

cryogenics

From cool-looking ice containers, to freezer racking compatible with most major brands, our innovative range of cryo products continues to expand.

Freezer Racks

Offering a variety of racks to safely store samples in your freezer

Choose from a selection of corrosion-resistant stainless steel racks that are compatible with most brands of upright and chest freezers.

- All have handles for easy use
- Racks store boxes up to a 135 x 135 mm footprint and accommodate 2" tall boxes (up to 5 cm) as well as 3" tall boxes (up to 7.6 cm) depending on the rack model

Chest Freezer Racks

Part no.	Config.	Box Size	L x W x H	Qty.
HS2862AB	7 x 1	2 in (5.3 cm)	5.5 x 5.5 x 15.3 in (14.1 x 14.1 x 39 cm)	1/ea
HS2862A	10 x 1	2 in (5.3 cm)	5.5 x 5.5 x 22 in (14.1 x 14.1 x 55.7 cm)	1/ea
HS2862AC	11 x 1	2 in (5.3 cm)	5.5 x 5.5 x 24.1 in (14.1 x 14.1 x 61.3 cm)	1/ea
HS2862AD	12 x 1	2 in (5.3 cm)	5.5 x 5.5 x 26.3 in (14.1 x 14.1 x 66.9 cm)	1/ea
HS2862AE	13 x 1	2 in (5.3 cm)	5.5 x 5.5 x 28.5 in (14.1 x 14.1 x 72.4 cm)	1/ea
HS2862DA	5 x 1	3 in (7.6 cm)	5.5 x 5.5 x 16 in (14.1 x 14.1 x 40.6 cm)	1/ea
HS2862D	8 x 1	3 in (7.6 cm)	5.5 x 5.5 x 25.5 in (14.1 x 14.1 x 64.9 cm)	1/ea
HS2862DB	9 x 1	3 in (7.6 cm)	5.5 x 5.5 x 29 in (14.1 x 14.1 x 73 cm)	1/ea

Upright Freezer Racks 2 inch Boxes

Part no.	Config.	Box Size	L x W x H	Qty.
HS2862HA	4 x 3	2 in (5.3 cm)	16.4 x 5.5 x 8.9 in (41.8 x 14.1 x 22.5 cm)	1/ea
HS2862H	4 x 4	2 in (5.3 cm)	21.8 x 5.5 x 8.9 in (55.6 x 14.1 x 22.5 cm)	1/ea
HS2862HD	4 x 5	2 in (5.3 cm)	27.3 x 5.5 x 8.9 in (69.5 x 14.1 x 22.5 cm)	1/ea
HS2862HE	5 x 3	2 in (5.3 cm)	16.4 x 5.5 x 11 in (41.8 x 14.1 x 28 cm)	1/ea
HS2862HB	5 x 4	2 in (5.3 cm)	22 x 5.5 x 11 in (55.6 x 14.1 x 28 cm)	1/ea

Upright Freezer Racks 3 inch Boxes

Part no.	Config.	Box Size	L x W x H	Qty.
HS2862KA	3 x 3	3 in (7.6 cm)	16.4 x 5.5 x 8.9 in (41.8 x 14.1 x 22.6 cm)	1/ea
HS2862K	3 x 4	3 in (7.6 cm)	22 x 5.5 x 8.9 in (55.6 x 14.1 x 22.6 cm)	1/ea
HS2862KB	4 x 4	3 in (7.6 cm)	22 x 5.5 x 12.8 in (55.6 x 14.1 x 32.5 cm)	1/ea
HS2862KC	5 x 4	3 in (7.6 cm)	22 x 5.5 x 16 in (55.6 x 14.1 x 40.6 cm)	1/ea



Handles for easy transfer

SEE PAGE 29-35
For storage boxes



Robust mounting brackets assure handles stay securely in place for easy use


4 x 3




4 x 4




4 x 5




5 x 3




5 x 4



3 x 3



3 x 4



Upright freezer rack configurations