





Product Features

Non-treated polystyrene dishes are transparent and DNase/RNase-Free. The bottom has a ridge allowing for easy handling and stacking. The plate has three vents for gas exchange.

- High clarity polystyrene
- Flat transparent surface
- PA/PE packaging
- Stackable for easy storage and handling
- Sterilized by E-Beam
- Non-pyrogenic



Cat. No.	Description	Height(mm)	Volume (ml)	Sterile	/Pack	/Case
752001	90mm	15	40	+	20	500
752011	90mm, 2 compartments	15	20	+	20	500
752021	90mm, 3 compartments	15	12	+	20	500
752031	90mm, 4 compartments	15	10	+	20	500
752002	90mm, Double bag package	15	40	+	20	500



Product Features

- Polypropylene material, hydrophilic
- Different colors for easier identification
- Flexible handles for easier collection and inoculation
- Sterilized by Ethylene Oxide
- Easy-tear individual packaging

Cat. No	Volume(ul)	Color	Length(mm)	/Pack	/Case
716001	-	White	225	400	4,000
717101	1	Blue	225	400	4,000
718201	10	Yellow	225	400	4,000