

NEW PRODUCT

NEW!



DIAMOND CULTURE TUBE CAPS

These popular caps fit plain end culture tubes and flasks with DeLong style necks or ports. The caps are molded from durable, autoclavable polypropylene (PP). They are available in a full range of five sizes and five colors for easy sample identification and organization. Recommended tube diameter is indicated on each cap.

- Polypropylene (PP)
- Autoclavable
- Ridges inside cap provide secure grip and gas exchange
- Available in five sizes and five colors
- Maintain sterility without cotton plugging
- Reduce sample evaporation compared to cotton plugging
- Provide same evaporation and sterility benefits as stainless steel closures



By popular demand!

Max Tube Diameter	Bag Qty.	Unit	Item#				
			Natural	Blue	Green	Red	Yellow
13mm	100	1000	118150C	118150B	118150G	118150R	118150Y
16mm	100	1000	118152C	118152B	118152G	118152R	118152Y
18mm	100	1000	118154C	118154B	118154G	118154R	118154Y
20mm	100	500	118156C	118156B	118156G	118156R	118156Y
25mm	100	500	118158C	118158B	118158G	118158R	118158Y



www.globescientific.com

Customer service, sales and technical support:

Phone 1 (800) 394-4562 • 1 (201) 599-1400

Fax 1 (201) 599-1406

E-mail mail@globescientific.com