

Benchmate C12V Centrifuge Features

User Friendly | High Speed | Durable



15000 rpm / 15596 x g Maximum Speed

Small Footprint

Saves valuable bench space

Lid Lock Safety

Lid opens automatically on run completion

Microprocessor Controlled

With Last Run Memory Feature

Imbalance Detection with Auto Cut-Off

Intuitive, Simple Interface

Quick and convenient settings

Brushless DC motor

Maintenance free, delivers better power/weight ratio

Large Digital Display

Back-Lit LCD for easy reading of multiple parameters

One Touch Short Spin Operation & RPM/RCF Conversion

USB port

For remote operation & data logger feature

User Programmable

User can set up to 99 programs of 4 lines, each as per run requirements



12 Place Microcentrifuge Tube Rotor
(Included)



Optional PCR Strip Rotor

2 X 8 X 0.2 mL PCR Tubes

Part # C12V-PCR-Rotor