

PD-Tip™ // Syringe Tips

For Repeating Pipettes

The Perfect Tips for Your Repeating Pipette.

BRAND® encoded PD-Tip™ // syringe tips are manufactured from high-quality virgin plastics for use with all standard repeating pipettes in a wide variety of applications.

- **Compatible with Most Pipettes:** BRAND® PD-Tip™ // Precision Dispenser Tips are compatible with repeating pipettes using standard tips. BRAND® PD-Tip™ // Precision Dispenser Tips are encoded for automatic recognition by the BRAND® HandyStep® touch repeating pipette, the Rainin® AutoRep™ E pipette, and the Gilson® Repetman®.
- **Reliable:** A precision seal between the piston and cylinder ensure smooth operation, increased accuracy, and increased precision. Extended piston design minimizes dead volume and maximizes fluid utilization.
- **Manufactured from High Quality Materials:** PD-Tip™ // syringe tips are made entirely from high-quality virgin polyethylene and polypropylene. (Piston of 0.1 mL tip is liquid crystal polymer-LCP).
- **BIO-CERT® PD-Tip™ // Precision Dispenser Tips** are sterile, individually wrapped, and certified free of DNA, RNase, endotoxins, and ATP to meet the most stringent sample integrity requirements.
- **Batch Certificate of Performance** which documents tip accuracy and precision testing is supplied with each pack of BRAND® PD-Tip™ // Precision Dispenser Tips.



PD-Tip™ // Precision Dispenser Tips BIO-CERT® sterile, individually wrapped

Capacity,	Pack of	Cat. No.
0.1 mL	100	705730
0.5 mL	100	705732
1 mL	100	705734
1.25 mL	100	705736
2.5 mL	100	705738
5 mL	100	705740
10 mL	100	705742
12.5 mL	100	705744
25 mL (one adapter included)	25	705746
50 mL (one adapter included)	25	705748
Adapter for 25 mL and 50 mL , BIO-CERT® sterile, PP, autoclavable		
25 mL and 50 mL Adapter	5	702399

*The Eppendorf® Repeater® Plus, Repeater® M4, Repeater® Pro, Stream, and XStream are designed to preclude their operation with any tips other than Eppendorf® Combitips® advanced.

PD-Tip™ // Precision Dispenser Tips, non-sterile

Capacity,	Pack of	Cat. No.
0.1 mL	100	705700
0.5 mL	100	705702
1 mL	100	705704
1.25 mL	100	705706
2.5 mL	100	705708
5 mL	100	705710
10 mL	100	705712
12.5 mL	100	705714
25 mL (one adapter included)	50	705716
50 mL (one adapter included)	25	705718
Multi-pack (20 each 0.5, 1, 1.25, 2.5, 5, 10 and 12.5 mL)		
Adapter for 25 mL and 50 mL , non-sterile, PP, autoclavable		
25 mL and 50 mL	10	702398



Lab Rats Trust BrandTech!

Toll Free: 888-522-2726 www.brandtech.com

HandyStep_BT_Fly-2021

BRANDTECH®
SCIENTIFIC, INC.