

Ready·Strain™ Cell Straining Kits



- Includes strainer, 50mL tube, & cap
- Opens directly under the hood reduces contamination risks
- Preassembled - ready to use
- Guaranteed sterile - gamma
- Packaged in easy-open “tear-strip” bags

ReadyStrain™ is the first product of its kind, providing all the necessary parts for filtration of cell suspensions and combining them into a **single, preassembled, ready to use, sterile packet**. These kits eliminate the contamination risks associated with opening multiple packages, and assembling separate pieces (i.e. tubes strainers and caps). The individual sterile packets mean researchers are no longer forced to expose entire bags of 50mL tubes and corresponding caps in order to perform one filtration. Instead they can open a single ReadyStrain packet directly under the hood.

Each kit contains a *sterile* SureStrain™ Strainer, pre-attached to a *sterile* 50mL tube, with a separate *sterile* screw-cap ready for assembly, processing and storage after the filtration is complete. Kits are supplied in easy open “tear-strip” bags and gamma sterilized.

Pipette.com

10360 Sorrento Valley Rd.
San Diego, CA 92121
Ph: 1800-243-3232 www.pipette.com

Effective Jan. 1, 2018. Subject to change without notice.

What's included:

1 x 50mL Centrifuge Tube w/ graduations

1 x Attached SureStrain Cell Strainer
(40, 70 or 100μm)

1 x Screw Cap for 50mL Centrifuge Tube

**Kits supplied sterile, preassembled,
in individual easy-open bags, 50/pk**



Also available:
PolyPestle™ Dual-sided Pestle



<u>Cat. No.</u>	<u>Description</u>	<u>Qty</u>
C5040	Cell straining kit w/ strainer, 50mL tube and screw-cap, 40μm, blue, sterile, individually wrapped	50
C5070	Cell straining kit w/ strainer, 50mL tube and screw-cap, 70μm, white, sterile, individually wrapped	50
C5100	Cell straining kit w/ strainer, 50mL tube and screw-cap, 100μm, yellow, sterile, individually wrapped	50
C4010	PolyPestle™ dual-sided pestle for strainers & 1.5mL microtubes, sterile, individually wrapped	50