

droplet™

Disposable sample organization and storage

Corrugated Sample Storage Tray

Corrugated storage trays are multi-purposed where they can be used with the included dividers to hold specimen containers or used for general lab product storage with the inserts removed.

- Made from corrugated polypropylene
- Moisture resistant, significantly reducing the chance of mold
- Lightweight and sturdy to hold up to 24 full specimen containers
- Stackable with other trays
- Removable dividers allow outer shell to be used for storage or transport of laboratory products
- Supplied flat-packed and includes: 10 trays and 10 dividers
- Disposable
- Store down to temperatures of -80°C


| Part No. | Wells | Well Dia. | L x W x H | Pkg. | Color |
|----------|-------|-----------|--|---------|--------------------------------|
| HS120375 | 24 | 58 mm | 16 x 10.24 x 3.5 in (40.6 x 26 x 8.9 cm) | 10/pack | White <input type="checkbox"/> |



Use for storage or transport of laboratory products



Removeable inserts



Trays are supplied "flat-packed"

Reducing the amount of space needed for storage



Droplet™ Sample Storage Tray

Sample organization for the busy lab


Lightweight and sturdy to hold up to 24 full specimen containers (specimen containers not included).

- Removable insert allows outer shell to be used for storage or transportation of laboratory products
- Stackable with other trays to support as many as four high
- Integrated plastic label holder enables quick and easy change identification (cards Included)
- Able to be cleaned with common disinfectants and germicides



Use for storage or transport of laboratory

| Part No. | Wells | Well Dia. | L x W x H | Pkg | Color |
|-----------------|-------|-----------|--|------|--------------------------------|
| HS120172 | 24 | 58 mm | 18.9 x 13 x 3.3 in (48 x 33.1 x 8.4 cm) | 1/ea | White <input type="checkbox"/> |



| Part No. | Array | Wells | Wells Dia. | L x W x H | Pkg. | Color |
|-----------------|-------|-------|------------|---|------|--|
| HS120091 | 6 x 4 | 24 | 58 mm | 15.47 x 10.31 x 1.97 in (39.3 x 26.2 x 5 cm) | 1/ea | Blue <input checked="" type="checkbox"/> |

Urine Container Rack

Durable rack for sample handling

- Ideal for urine containers up to 58 mm dia.
- Slotted bottom wire partitions
- Mesh bottom support
- Autoclavable