

## NEW! 15mL Polystyrene Centrifuge Tubes, Sure Cap

- Polystyrene tube is optically clear for easy viewing of sample
- Sure Cap with compressible liner creates a stable seal
- Includes a molded 15mL graduation line and 0.5 and 1mL graduations molded into the conical bottom allowing for small volume quantification
- Frosted marking area
- Max Force (RCF) is 3,000g



CELLTREAT Part No.	Product Description	Pack Description	Qty/Case
229491	15mL Polystyrene Centrifuge Tube, Sure Cap	25/Bag	500

## NEW! 15 & 50mL Centrifuge Tubes, Sure Cap

- Sure Cap with compressible liner creates a stable seal
- Graduations are printed in bold easy-to-read black ink and are accurate to  $\pm 2\%$
- Large white marking area provides excellent contrast for identification
- Medical grade virgin autoclavable polypropylene & have smooth, flat surface polyethylene (HDPE) caps
- Temperature range of  $-80^{\circ}\text{C}$  to  $120^{\circ}\text{C}$
- Sterile (gamma irradiated) in zip-closure bags



CELLTREAT Part No.	Product Description	Pack Description	Qty/Case
229401	15mL Centrifuge Tube, Sure Cap	25/Bag	500
229402	50mL Centrifuge Tube, Sure Cap	25/Bag	300