

SureStand™

MULTI-CHANNEL CAPABLE PIPETTE RACK



- **Universal - Fits all popular pipettes**
- **5mm acrylic - 40% thicker**
- **Non-skid rubber feet**
- **Multi-channel “window” slots**

P4403



3-Place Rack
up to one multi-channel

P4405




5-Place Rack
up to two multi-channels

P4406

SureStand™ Pipette racks are unique in their ability to accommodate multi-channel pipettes without sacrificing any of the stability of a traditional acrylic rack. The patent pending “windows” allow the user to insert the multi-channel head longways through the rack, while keeping the body of the pipette stable by securing the pipette at two points; under the finger hook, and near the bottom of the handle.

SureStand is available in 3 versions; a 3-place rack that accepts one multi-channel, a 5-place rack that accepts two multi-channels, and a 6-place rack that accepts four multichannels. Each SureStand is also outfitted with four non-skid rubber “feet” for extra stability. SureStand is guaranteed to be the sturdiest and most versatile pipette rack in your lab



SureStand’s “window” slots are the first of their kind, allowing the user to store single channel or multi-channel pipettes easily and conveniently.

MTC Bio

2600A Main Street Extension
Sayreville, NJ 08872
Ph: 908-279-0280 www.mtcbio.com

Ordering Information

<u>Cat. No.</u>	<u>Description</u>	<u>Price</u>
P4403	SureStand Pipette Rack for 3 pipettes (up to one multi-channel)	\$45.00
P4405	SureStand Pipette Rack for 5 pipettes (up to two multi-channels)	\$69.00
P4406	SureStand Pipette Rack for 6 pipettes (up to four multi-channels)	\$77.00