

DEEP WELL PLATES



T110-2 & -3

BioBlock™ Deep Well Plates with 600 µl 8-Tube Strips



Made of polypropylene


These plates feature 600 µl wells for smaller volume applications. They include a T110-5 deep well plate along with twelve T110-15 eight-tube strips (see below). In the T110-2 model, tube strips are removable and can also be ordered separately. In the T110-3 model, tube strips are welded by ultrasound.

DMSO resistant. Packed in bags of 4 plates.



Cat. #	Description	Qty/Pk	Qty/Cs
T110-2	Plate, 600 µl (removable tube strips)	4	24
T110-3	Plate, 600 µl (fixed tube strips)	4	24



96 Wells with Conical Bottom



600 µl
WORKING
VOLUME



T321-1 & -2

Domed and Flat Cap Strips

Made of polypropylene

For a perfect seal, 8-cap strips are available.

Cat. #	Description	Color	Qty/Cs
T321-1N	Domed cap strip	Natural	125
T321-2N	Flat cap strip	Natural	125



T110-15

Strip of 8 Tubes 600 µl

Made of polypropylene

Cat. #	Description	Qty/Cs
T110-15	Strip of 8 tubes, 600 µl	125