

MTC Bio™

DISPOSABLE TRANSFER PIPETTES



- **Transparent, nontoxic & inert**
- **Precise, uniform drop size**
- **Shatter-proof & unbreakable**
- **Can be used with LN₂ (-196°C)**
- **Heat sealable**

These unbreakable, one-piece pipettes eliminate contaminated rubber bulbs, broken glass and exposure to infectious materials. They are made of ultra-clear, low density polyethylene (LDPE) and are ideal for microbiology, blood banking, hematology, immunology, and urinalysis. Sterile pipettes are excellent for tissue culture applications. Their low-affinity surface reduces sample, cell or protein loss.

Whether you make your selection based on length, bulb draw, pipette capacity or number of drops per ml, one of these ten designs is likely to be suitable for your application. Some have ultra-fine extended tips (P4113, P4117, P4121, P4122, P4133) that allow for extremely precise dispensing and access into the narrowest vessels. They can also be heat sealed for airtight sample storage. Their uniform drop size can be used for reproducible spotting of gels & TLC plates, and for dispensing into test plates, ID strips, PCR tubes or 96 well plates.



2600A Main Street Extension
Sayreville, NJ 08872
Ph: 908-279-0280 www.mtcbio.com

Ordering Information

Cat. No.	Description	Packaging	Qty/pk	Price
P4111-00		Nonsterile, bulk	250	\$12.50
P4111-11		Sterile, Individ. Wrap	500	\$49.50
P4111-14		Sterile, 20/pk	500	\$39.00
P4114-00		Nonsterile, bulk	250	\$12.50
P4114-11		Sterile, Individ. Wrap	500	\$49.50
P4114-14		Sterile, 20/pk	500	\$39.00
P4117-00		Nonsterile, bulk	250	\$12.50
P4117-11		Sterile, Individ. Wrap	500	\$49.50
P4117-14		Sterile, 20/pk	1500	\$39.00
P4122-01		Nonsterile, bulk	500	\$45.00
P4122-11		Sterile, Individ. Wrap	500	\$64.00
P4131-00		Nonsterile, bulk	250	\$12.50
P4131-11		Sterile, Individ. Wrap	500	\$49.50
P4131-14		Sterile, 20/pk	500	\$39.00

Cat. No.	Description	Packaging	Qty/pk	Price
P4113-00		Nonsterile, bulk	250	\$12.50
P4113-11		Sterile, Individ. Wrap	500	\$49.50
P4115-00		Nonsterile, bulk	250	\$12.50
P4115-11		Sterile, Individ. Wrap	500	\$49.50
P4115-14		Sterile, 20/pk	500	\$39.00
P4121-01		Nonsterile, bulk	500	\$45.00
P4121-11		Sterile, Individ. Wrap	500	\$64.00
P4123-01		Nonsterile, bulk	400	\$45.00
P4123-11		Sterile, Individ. Wrap	500	\$64.00
P4123-14		Sterile, 20/pk	500	\$62.00
P4133-00		Nonsterile, bulk	250	\$12.50
P4133-11		Sterile, Individ. Wrap	500	\$49.50
P4133-14		Sterile, 20/pk	500	\$39.00