

# AccuPet Jr.



one precise task - two quality solutions



# AccuPet Pro

# AccuPet Jr.

AccuPet Jr. is the first two step plunger dispensing and smart tip cone design color coded minipipette that provides the same level of reliability as standard pipettes Super light and fully autoclavable pipette provides superb accuracy at the fraction of the cost of average fixed volume pipette . Perfect for diagnostics kits, advertisement as well as many other measuring applications. Pipette may be supplied with many customizable options.



# AccuPet Pro Fixed

Designed to minimize the amount of hand stress over time, lightweight pipette feature rubberized soft surfaces that helps to make pipetting more comfortable than never before. Accuracy and precision improves while reducing the risk of cumulative trauma disorders. Pipette feature universal tip fitting shaft, light ejection mechanism and come with calibration certificate and three year warranty.



Model	Volume (µl)	Color	Accuracy (%)	Precision (%)	List Price (\$)
AJ-5	5.0	Red	1.6%	1.0%	\$60 (set of 5)
AJ-10	10.0	Pink	1.6%	1.0%	\$60 (set of 5)
AJ-20	20.0	Green	1.2%	0.8%	\$60 (set of 5)
AJ-25	25.0	Green	1.2%	0.8%	\$60 (set of 5)
AJ-50	50.0	Purple	1.2%	0.8%	\$60 (set of 5)
AJ-100	100.0	Red	1.2%	0.8%	\$60 (set of 5)
AF-5	5.0	Red	0.8%	0.4 %	\$129.00
AF-10	10.0	Red	0.6%	0.3%	\$129.00
AF-20	20.0	Yellow	0.6%	0.3%	\$129.00
AF-25	25.0	Yellow	0.6%	0.3%	\$129.00
AF-50	50	Blue	0.6%	0.2%	\$129.00
AF-100	100	Purple	0.5%	0.2%	\$129.00
AF-200	200.0	Purple	0.5%	0.2%	\$129.00
AF-250	250.0	Purple	0.4%	0.2%	\$129.00
AF-500	500.0	Navy	0.4%	0.2%	\$129.00
AF-1000	1000.0	Navy	0.3%	0.08%	\$129.00
AF-2000	2000.0	Teal	0.3%	0.06%	\$129.00
AF-5000	5000.0	Teal	0.3%	0.04%	\$129.00
AF-CST / AJ-CST	Custom volume sizes are available upon request				



**Advanced Pipetting Systems**  
 13354 Russet Leaf Lane  
 San Diego • CA 92129  
 800-455-1779 (p) • 858-536-3342 (f)  
 www.pipetpro.com • info@pipetpro.com

Distributed by: