

NEW!



WIZARD MINI-FUGE

Compact and Space-Saving Personal Centrifuge!



WZ-MF6000

- Sleek and Compact mini centrifuge is ideal for microfiltration, HPLC samples, DNA, RNA, and Protein applications.
- Includes adapters to accommodate microtube and PCR tube styles (2.0, 1.5, 0.5, 0.2mL) for customized use.
- Lid functionally designed to fit in the palm of the hand for easy repeat opening and quick tube extraction.
- Tool-free snapspin™ rotors and quick-spin options make it fast and easy to change out rotors with a push of a button..
- Weighted rubberized base eliminates movement and dampens vibration and noise.
- Includes a tube rack and storage case for rotors and adapters.
- The WIZARD Mini-Fuge is a perfect compliment to the Vortex-Genie 2 Family.



Accessories Included

Catalog No.

Description

WZ-MF6000

Wizard Mini-Fuge

Specifications

Supplied with

6 1.5/2.0ml tube rotor, 16 place 0.2ml tube rotor (singles, or 2x8 strip), 6 tube adapters 0.2ml, 6 tube adapters 0.5ml, 4 o-rings, storage case, tube rack

Maximum Speed

6000 RPM/2000 x g

Dimensions (D X W X H)

153 x 128 x 104mm (6 x 5 x 4in)

Electrical

100-240V, 50/60 Hz 0.6A Approx. 7.3W (supplied with 4 interchangeable power plugs)



80 Orville Drive, Suite 102
Bohemia, NY 11716 U.S.A.
Phone: (631) 567-4700 • Fax: (631) 567-5896
Toll Free: 888-850-6208
Customer Service: custsvc@scientificindustries.com
scientificindustries.com